

# ENDO

## ENDO Baby Incubator EI.BI



### Spesifikasi

<b>Power supply</b>	AC110/220V 60/50Hz	<b>Power input</b>	650 VA
<b>Kondisi Pengoperasian</b>			
<b>Suhu</b>	20~30°	<b>Ukuran kasur</b>	65cm (L) x 37cm (W)
<b>Kelembaban relative</b>	30~75%	<b>Kapasitas Palung</b>	1000 mL
<b>Kecepatan aliran udara</b>	<0.3 m/detik	<b>Penyaring udara</b>	0.5 µm
<b>Rentang kontrol suhu udara</b>	25°C ~ 37°C 37.1°C ~ 38°C	<b>Dimensi</b> Setiap unit dikemas dalam satu kotak polywood (unit kontrol dan alas)	126 x 71 x 85 cm
<b>Fluktuasi suhu</b>	±5°C		
<b>Keseragaman suhu kasur</b>	≤0.8°C	<b>Berat</b>	Net 80Kg Gross 100Kg
<b>Tingkat kebisingan internal</b>	<55dB (A)		
<b>Sudut kemiringan tempat tidur bayi</b>	±12°		