

ENDO

ENDO Chemistry Analyzer EI.CA2



Spesifikasi

Wavelength (Panjang gelombang)	340, 380, 405, 505, 546, 578, 620 (nm) dan satu posisi bebas
Light source	Halogen lamp, 6V 10W
Tampilan	5 inch touchscreen
Memory	15.000 hasil penyimpanan
Volume aspirasi	100~5000 μ L
Metode test	End Point, Fixed time, Kinetics, SUB, DIV, ABS
Reagen	Buka system reagen
Rentang absorbansi	0~3.000 Abs
Variasi absorbansi (340 nm)	\leq 0.005 Abs
Ukuran	355 x 330 x 175 mm
Power supply	AC 100-240V 50/60HZ
Power	\leq 180 VA
Berat bersih	5.6 kg