


ENDO



ENDO EI.AM7 **Anesthesia Machine**

<https://endo.id>

ENDO EI.AM7

Anesthesia Machine

SPEKIFIKASI TEKNIS

CATU DAYA AC EKSTERNAL

Tegangan Input	100~240 VAC
Frekuensi Input	50/60Hz
Arus input maksimum	6.3A

BATERAI INTERNAL


Jumlah Baterai	1
Jenis Baterai	Li-ion Batter
Tegangan Nominal Baterai	14.8V DC
Kapasitas Baterai	5200mAh

SPEKIFIKASI PERANGKAT KERAS


Layar	TFT LCD berwarna
Antarmuka	<ul style="list-style-type: none">•Daya : 1 input AC + 3 output daya bantu•Jaringan : 1 port Ethernet•Monitor : VGA 15-pin•RS232 : 1 port•Grounding : 1 terminal pentanahan equipotensial
Mobilitas	Roda: 4 buah (diameter 12,5 cm), semua dilengkapi rem

ENDO

 Jalan Raya Menganti 14 Kedurus, Surabaya 60223
Jawa Timur, Indonesia

 +62-31-767-3636

 info@endo.id

 +62-31-767-3737

 <https://endo.id>

 0800 . 177 . ENDO (3636)