


# ENDO Portable Endoscopy Camera System EC-1



CE FDA

### Product Parameter

- ◆ Camera: 1/3" CCD waterproof camera
- ◆ Scanning standard: DVI output 600 lines
- ◆ Output: Video
- ◆ Light source: LED lamp
- ◆ Light source color temperature: 6000K
- ◆ Storage: SD card dynamic/ static storage
- ◆ Display screen: 15" TFT color screen
- ◆ Screen resolution: 800 x 600



### ◆ EASY AND PRACTICAL OPERATION

- Simple panel adjustment can get a nice image effects
- White balance support manual or automatic adjustment

### ◆ APPLIES TO ALL CLINICAL DEPARTMENTS

- General Surgery
  - Gynecology
  - ENT
  - Gastroenterology
  - Respiratory Medicine
  - Urology
  - Orthopedics
  - Thoracic Surgery
- All departments using the endoscope

## PT. ENDO INDONESIA

Jl. Raya Menganti No. 14, Wiyung - Surabaya 60223

Tel : (031) 767 3636

Email : [info@endo.id](mailto:info@endo.id)

Fax : (031) 767 3737

URL : <https://endo.id>