

# ENDO

## ENDO Electrosurgical EI.ESU



### Spesifikasi

Nilai tegangan	AC120V / AC230V
Frekuensi	50Hz / 60Hz
Konsumsi daya	1200VA + 10%
Frekuensi Carier	400kHz, 500kHz
Frekuensi pengulangan	33kHz
Ukuran	330 x 455 x 150mm
Berat	13Kg

### Spesifikasi Output

Pure cut	400W / 300Ω
Blend 1	250W / 300Ω
Blend 2	200W / 300Ω
Blend 3	150W / 300Ω
Contact Coag.	120W / 300Ω
Spray Coag.	100W / 300Ω
Bipolar Coag.	100W / 100Ω
Bipolar Cut	100W / 200Ω