

ENDO

ENDO Electrosurgical EI.ESU.10



Spesifikasi

| | |
|------------------|-----------------------|
| Nilai tegangan | AC120V ~ AC230V |
| Frekuensi | 50Hz / 60Hz |
| Konsumsi daya | 1100VA+10% |
| Ukuran | 330x455x150mm |
| Berat | 13Kg |
| Suhu operasi | 10°C hingga 40°C |
| Suhu penyimpanan | -10°C hingga 60°C |
| Kelembaban | 20% hingga 95% RH |
| Tekanan | 700~1060mbar |
| Siklus operasi | 10 dtk ON 30 dtk diam |

Output Power

| | | |
|----------|---------------|-------------|
| Pure cut | Output power | 300W / 300Ω |
| | Carrier Freq | 400kHz |
| | Crest Factors | 1.6 |
| | Duty Rate | 100% |
| | Repeat Freq | Continuous |
| Blend 1 | Output power | 230W / 300Ω |
| | Carrier Freq | 400kHz |
| | Crest Factors | 2.0 |
| | Duty Rate | 80% |
| | Repeat Freq | 33kHz |
| Blend 2 | Output power | 180W / 300Ω |
| | Carrier Freq | 400kHz |
| | Crest Factors | 2.2 |
| | Duty Rate | 60% |
| | Repeat Freq | 33kHz |
| Blend 3 | Output power | 120W / 300Ω |
| | Carrier Freq | 400kHz |
| | Crest Factors | 2.4 |
| | Duty Rate | 50% |
| | Repeat Freq | 33kHz |

Output Power

| | | |
|---------------------|---------------|-------------|
| Contact Coagulation | Output power | 100W / 300Ω |
| | Carrier Freq | 400kHz |
| | Crest Factors | 3.3~1.5 |
| | Duty Rate | 10% |
| | Repeat Freq | 33~153kHz |
| Spray Coagulation | Output power | 80W / 300Ω |
| | Carrier Freq | 400kHz |
| | Crest Factors | 4.5 |
| | Duty Rate | 8.0% |
| | Repeat Freq | 33kHz |
| Bipolar Coagulation | Output power | 80W / 100Ω |
| | Carrier Freq | 500kHz |
| | Crest Factors | 1.5 |
| | Duty Rate | 100% |
| | Repeat Freq | Continuous |

Optional (Upgrade)

| | | |
|-------------|---------------|------------|
| Bipolar Cut | Output power | 80W / 200Ω |
| | Carrier Freq | 500kHz |
| | Crest Factors | 1.5 |
| | Duty Rate | 100% |
| | Repeat Freq | Continuous |