

ENDO

ENDO Electrosurgical EI.ESU2



Spesifikasi

Nilai tegangan	AC120 / AC230V
Frekuensi	50Hz / 60Hz
Konsumsi daya	1200VA+10%
Frekuensi carier	400kHz, 500kHz
Frekuensi pengulangan	33kHz
Ukuran	390x410x160mm
Berat	12.5Kg
Suhu operasi	15°C hingga 30°C
Suhu penyimpanan	-10°C hingga 60°C
Kelembaban	20% hingga 95% RH,
Ketinggian operasi	700mbar~1060mbar
Siklus operasi	10 detik ON 30 detik Diam

Spesifikasi Output

Cut	Pure	400W / 200Ω
	Blend 1	250W / 200Ω
	Blend 2	200W / 200Ω
	Blend 3	150W / 200Ω
Coagulator	Contact	120W / 200Ω
	Spray	100W / 300Ω
	Argon	100W / 300Ω
Bipolar Coag	Standard	100W / 100Ω
	Auto Start	100W / 100Ω
	Forced	80W / 100Ω
Bipolar Cut	Cut	120W / 200Ω
	Blend	100W / 200Ω
Option 1		
Endo Cut	Papilotomy	150W / 100Ω
	Polypectomy	150W / 100Ω
Coagulation	Soft	80W / 100Ω