

ENDO

ENDO Electrosurgical EI.ESU.9



Spesifikasi

Nilai tegangan	AC120V ~ AC230V
Frekuensi	50Hz / 60Hz
Konsumsi daya	550VA+10%
Ukuran	330x330x135mm
Berat	8.5Kg
Suhu operasi	10°C hingga 40°C
Suhu penyimpanan	-10°C hingga 60°C
Kelembaban	20% hingga 95% RH
Tekanan	700~1060mbar
Siklus operasi	10 dtk ON 30 dtk diam

Output Power

Pure cut	Output power	150W / 200Ω
	Carrier Freq	400kHz
	Crest Factors	1.6
	Duty Rate	100%
	Repeat Freq	Continuous
Blend 1	Output power	120W / 200Ω
	Carrier Freq	400kHz
	Crest Factors	1.9
	Duty Rate	67%
	Repeat Freq	33kHz
Blend 2	Output power	100W / 200Ω
	Carrier Freq	400kHz
	Crest Factors	2.3
	Duty Rate	42%
	Repeat Freq	33kHz
Contact Coagulation	Output power	80W / 200Ω
	Carrier Freq	400kHz
	Crest Factors	3.0
	Duty Rate	25%
	Repeat Freq	33kHz
Spray Coagulation	Output power	60W / 500Ω
	Carrier Freq	400kHz
	Crest Factors	4.1
	Duty Rate	17%
	Repeat Freq	33kHz
Bipolar Coagulation	Output power	80W / 100Ω
	Carrier Freq	500kHz
	Crest Factors	1.6
	Duty Rate	100%
	Repeat Freq	Continuous