

ENDO

ENDO Infusion Pump EI.IP3



Spesifikasi

Dimensi	415 mm x 245 mm x 215 mm
Berat	1.7 kg
Layar	Display 5.0" ; Touchscreen
Resolusi	480*320
Kecerahan	Dapat disesuaikan, tingkat 1-10
Drop sensor	Standar

Spesifikasi

Spesifikasi Infus	Mode infus	Mode Rate, Time (Waktu), VTBI, BW, Intermittent, Loading Dose, Sequence, Trapezia, Drip
	Laju aliran (Flow rate)	0.10-1600 ml/jam
	Kenaikan (Increment)	0.10-99.99 ml/jam : 0.01 100.0-999.9 ml/jam : 0.1 1000-1600 ml/jam : 1
	Volume preset (VTBI)	0-9999.99 ml
	Preset time	00:00:01-99:59:59 (jam : menit : detik)
	KVO	0.10-5.00 ml/jam Kenaikan : 0.01 ml/jam Default 1.00 ml/jam
	Purge	0.10-1600 ml/jam
	Akurasi	± 5%
Bolus	Bolus rater	Bolus Auto/manual 0.10-1600 ml/jam
	Preset volume bolus	0.10-50.00 ml Kenaikan : 0.01 ml
Sistem keamanan	Tipe alarm	Alarm visual dan akustik
	Volume	Tingkat 1-8
	Tipe log peristiwa	Operasi, Alarm, Infus
	Nomor maks	4000
	Preset nomor obat	80