

ENDO

ENDO Nebulizer EI.Neb02



Spesifikasi

Sumber daya	AC230V, 50Hz
Nilai daya	≤140VA
Aliran udara	10-12L/menit
Noise	≤65dB (A)
Tingkat semprotan	≥0.15ml/menit (NaCl 0.9%)
MMAD	partikel teratomisasi berdiameter 1μm-5μm menyumbang > 70%
Kapasitas wadah	2-8ml
Residu obat	≤1ml
Ukuran	266mm(L)x183mm(W)x128mm(H)
Berat	1.9Kg