

# HANDPIECE AKD

## ENDO EI.PT.BD

Bone Drill

KEMENKES RI AKD 21603220003



### Features

- ENDO Bone Drill EI.PT.BD is an instrument used for hole drilling applications for plates and screws during trauma surgery. This instrument has small volume, light weight, no power cord, easy to carry and operate, thereby reducing operating time.

### Specification

- Voltage : 14 V
- Power : 80W
- Temperature : < 25°C
- Output power : AC 110-220V, 50/60Hz 35va
- Output Voltage : 14V, 700Ma
- Charging time : < 3 hour
- Input Power : 30VA
- Battery : 14.4V
- Speed : 1100 r.p.m

### Bone Drill Set



Bone Drill  
Quantity : 1



Battery  
Quantity : 2



Charger  
Quantity : 1



Key  
Quantity : 1



Sterilized case pckaging  
Quantity : 1