

# HANDPIECE AKD

## EI.PT.CD V2

Cannulated Drill

KEMENKES RI AKD 21603220504



### Features

- The ENDO Cannulated Drill (EI.PT.CD.v2) is a medical electric drill consisting of a handpiece, charger, battery, and accessories.
- The housing/package of the unit is completely enclosed and can be sterilized in an autoclave.
- With its ergonomic design, the entire unit is autoclavable, low noise, high speed, alignment, and long service life.

### Specification

- Voltage : 14 V
- Output Power :  $\geq 100$  W
- Temperature :  $\leq 25^{\circ}\text{C}$
- Output Power : AC 110-220 V, 50/60 Hz
- Output Voltage : 14.4 V, 700 mA
- Charging Time :  $\leq 3$  hours
- Input Power : 30 VA
- Battery : 14.4 V

### Cannulated Drill Set



Cannulated Drill  
Quantity : 1



Battery  
Quantity : 2



Charger  
Quantity : 1



Transfer Kit  
Quantity : 1



Key  
Quantity : 1



Sterilized case pckaging  
Quantity : 1