

ENDO EI.PT.MF

Multifunction Handpiece

KEMENKES RI AKD 21603323872



Features

- Orthopedic handpiece/power tool with a multifunction drill that can be used for trauma surgery (plate and screw implant placement), Kirschner wire placement (K-wire pinning), and acetabular reamer.
- It has four interchangeable attachment tips: AO Drill, Bone Drill, K-Wire Driver, and Acetabular Reamer.

Specification

- Voltage : 9.6V
- Output power : $\geq 100W$
- Noise : ≤ 55 dB (A)
- Motor Speed : ≥ 1200 RPM
- Output power : AC 110-240V, 50/60Hz
- Output voltage : 9.6V, 700mA
- Charging time : ≤ 1.5 hours
- Input power : 30VA
- Battery : Ni-MH 2/3A 1800mAh 9.6V

Cannulated Drill Set



Cannulated Drill
Quantity : 1



Quick Coupling
Quantity : 1



Acetabulum Reaming
Drill Attachment
Quantity : 1



AO Drill Attachment
Quantity : 1



Bone Drill Attachment
Quantity : 1



K-Wire Driver Attachment
Quantity : 1



Battery
Quantity : 2



Charger
Quantity : 1



Battery Case
Quantity : 2



Aseptic Transfer Kit
Quantity : 2



Key
Quantity : 1



Sterilized case packaging
Quantity : 1