

# ENDO

## ENDO Infrared Thermometer (Non-Touch) EI.TM.3



### Spesifikasi

Display	34 x 21mm LCD display
Akurasi Suhu telinga, suhu dahi	(35.0°~42.0°) ±0.2°; (95.0°~107.6°) ±0.4° (32.0°~34.9°) ±0.3°; (89.6°~94.9°) ±0.5° (42.1°~42.9°)±0.3°; (107.7°~109.2°) ±0.5°
Akurasi Suhu objek, suhu lingkungan	(0.0°~100.0°) ±1.0°; (32.0°~212.0°) ±1.8°

## Spesifikasi

Lingkup pengukuran Suhu dahi	32°~42.9° (89.6°~109.2°)
Lingkup pengukuran Suhu objek, suhu lingkungan	0°~100.0° (32.0°~212.0°)
Resolusi	0.1° (0.1°)
Pengukuran suhu	Infrared scan
Memory record	20 set
Power dan tegangan	DC 3V (2pcs 1.5V AAA battery)
Berat	Sekitar 57g (tidak termasuk baterai)
Dimensi	150 x 37 x 50mm