

ENDO

ENDO Digital Color USG EI.USG4D



TERANG

Pencahayaan yang baik saat pemindaian sehingga mata tidak mudah lelah

SUNYI

Probe bebas getaran dan bebas noise untuk merawat ibu dan janin dengan baik

PINTAR

Teknologi terbaru, i-image, SRA, dll. yang memberikan kualitas gambar super

CEPAT

Kecepatan volume 4D menangkap pergerakan janin secara real time

Spesifikasi

Mode Pencitraan

B, 2B, 4B, B/M, M

CFM

PW Mode

Power Doppler / Directional PD

Duplex

Trapezoidal

Real-time 4D (opsional)

Chroma B / PW

Probes

3.5 MHz Convex Probes

7.5 MHz Linear Probes

6.0 MHz Transvaginal Probe

4.0 MHz 4D Volume Probe

7.0 MHz Transvaginal Probe (180°)

3.0 MHz Micro-convex Probe Adult

5.0 MHz Pediatric Probe

Wideband, Multi-frequency

Opsional

7.5 MHz Linear Probe

6.0 MHz Transvaginal Probe

7.0 MHz Transvaginal Probe (180°)

3.0 MHz Micro-convex Probe (adult)

5.0 MHz Pediatric Probe

Video printer Sony UP897MD

PC printer (semua USB PC printer)

DICOM 3.0

Biopsy kit : convex, linear, TV probe

Foot-switch

Teknologi Pencitraan

THI

Speckle Reduction Algorithm

Compound Technology

i-Image : gambar optimal

Bidang Aplikasi

Obgyn

Vascular

Urologi

Basic Cardiac

Pediatric

Urologi

Muscoloskeletal

etc.

Bidang Aplikasi

Abdominal

Obstetrics

Gynecology

Vascular dan bagian kecil

Standar Konfigurasi

Unit utama, 19" LCD, 4 probe konektor

Hard disk (32GB), DVD-RW, 6 USB ports



3.5MHz Convex
D3C60L



7.5MHz Linear
D7L40L



6.0MHz Transvaginal
D6C12L



4.0MHz 4D Volume
V4C40L



7.0MHz Transvaginal
D7C10L(180°)



3.0MHz Micro-Convex
D3C20L



5.0MHz Pediatric
D5C20L

