



PM-1
Patient Monitor

About ENDO

ENDO INDONESIA was founded in 2006 as a company engaged in importation and distribution of Medical Equipment. Today our products include Electromedical Equipment, Reagents, Reusable and Disposable Medical Devices, among many others and accessories. And now we had grown into the manufacturing field of Medical Recording Chart Paper.

ENDO is a company that fulfills legal requirements and receives approval from company's legal entity of Indonesia, with number C-05336 HT.01.01.TH.2006. Supported by experienced technicians and proper equipments, we also provide repairs and calibration for electromedical equipment. ENDO exist to serve your various needs of Medical Equipments and highly concerned about customer satisfaction to promote long term working relationship.



PT. ENDO INDONESIA
Jl. Raya Menganti No. 14, Wiyung - Surabaya 60223
Tel : (031) 767 3636 | Email : info@endo.id
Fax : (031) 767 3737 | URL : <https://endo.id>
Toll Free : 0800 177 ENDO (3636)



ENDO PM-1

Patient Monitor



ENDO PM-1 Series

Reliable Patient Monitor



The ENDO PM-1 series put reliable monitoring technologies at your finger tips. With their accurate readings and optimized performance, the ENDO PM-1 series can meet basic monitoring needs, bringing easy access to high-quality healthcare for everyone.



ENDO PM-1 (10.4")

End tidal Carbon Dioxide for Intubated/Non-Intubated Patients



- ENDO G2 CO₂ (Sidestream)
 - Superior water trap design for accurate monitoring
 - iCARB™ algorithm with intelligent CO₂ pseudo wave identification technology
 - Multiple sampling accessories as options for adult, child and neonate patients
- Respironics CO₂ (Mainstream/Sidestream)
 - Plug & play module design
 - Dehumidification tube instead of water trap
 - Low sampling rate of 50ml/min suitable for all types of patients

Display Modes



3/5-lead ECG Analysis



OxyCRG for Neonate



Large Font



Trend Screen

Built-in Temporary Memory

500	96 h	60	720 s
NIBP Measurement Review	Trend Review	Alarm Review	Frozen Waveform

Seamless Connectivity

The Wi-Fi connectivity of the ENDO PM-1 series enables healthcare providers to monitor their patients' health status from almost anywhere. Connecting it with ENDO MFM-CMS software, you may log on from any PC / Tablet to and check the status of your patients.



- High resolution color TFT-LCD screen
- Pacemaker detection
- Defibrillation and electrosurgical interference protection
- Pitch tone (Pulse-tone modulation)
- Unique iSEAP™ algorithm which is specially optimized for arrhythmia patients
- iMAT™ SpO₂ algorithm with outstanding motion resistance and low perfusion resistance performance
- SP10 Validated NIBP
- Built-in rechargeable Li-Ion battery
- Nurse call
- VGA output (optional)
- Suitable for adult, pediatric and neonatal patients

Standard Parameters: 3/5-lead ECG, RESP, ENDO SpO₂, NIBP, PR, 2-TEMP
 Optional: 2-IBP, Respironics CO₂, ENDO G2 CO₂, Thermal Recorder, WLAN Accessory Kit



Rolling Stand

Wall Mount