


ENDO



ENDO EI.DFU

Disposable Flexible Video Ureteroscope

<https://endo.id>

ENDO EI.DFU

Disposable Flexible Video Ureteroscope

Single use flexible video ureterorenoscope delivers high quality image for visual inspection with a clean and safe clinical flow. It's affordable without reprocessing and maintenance, and it benefits both patient and hospital by effectively avoiding cross contamination.

Product Features

- 4-way bending facilitates multidimensional observation
- Excellent image quality with highly integrated CMOS and LED
- 3.6Fr working channel ensures good perfusion flow and reduces the risk of surgical infection
- Angle-lock function by resistance force gives precise angle control and reduces operating fatigue



Compound Bending
Up285°/Down285°



Slim Diameter



① **Capture/Recording**

② **Freeze**



Unique Waterproof Design



Left/Right Rotation

Specification

Model	Insertion Tube	Working Channel	Working Length	Field of View	Depth of Field	Bending Up/Down	Rotation (Left & Right)	Bending Logical type
URD-1S	7.5Fr	3.6Fr	670mm	90°	3-100mm	285°	120°	

Image Processor



*High Integration,
Portable, Touch-screen*




*1280*800
High Resolution*




*170 Wide-angle
Observation*



*All in one
Highly integrated design*



- Lag-free, ultra-sensitive touch operation
- Allow operating with gloves



Detachable rechargeable battery
supporting 4 hours working



Detachable SD Card



HDMI + BNC (CVBS) output

ENDO

📍 Jalan Raya Menganti 14 Kedurus, Surabaya 60223
Jawa Timur, Indonesia

☎ 62-31-7673636

☎ 62-31-7673737

☎ 0800 . 177 . ENDO (3636)

✉ info@endo.id

🌐 <https://endo.id>