



ENDO Sonotrax FD-1

Ultrasonic Pocket Fetal Doppler

About ENDO

ENDO INDONESIA was founded in 2006 as a company engaged in importation and distribution of Medical Equipment. Today our products include Electromedical Equipment, Reagents, Reusable and Disposable Medical Devices, among many others and accessories. And now we had grown into the manufacturing field of Medical Recording Chart Paper.

ENDO is a company that fulfills legal requirements and receives approval from company's legal entity of Indonesia, with number C-05336 HT.01.01.TH.2006. Supported by experienced technicians and proper equipments, we also provide repairs and calibration for electromedical equipment. ENDO exist to serve your various needs of Medical Equipments and highly concerned about customer satisfaction to promote long term working relationship.



PT. ENDO INDONESIA
 Jl. Raya Menganti No. 14, Wiyung - Surabaya 60223
 Tel : (031) 767 3636 Email : info@endo.id
 Fax : (031) 767 3737 URL : <https://endo.id>
 Toll Free : 0800 177 ENDO (3636)



ENDO Sonotrax FD-1

Ultrasonic Pocket Fetal Doppler



ENDO Sonotrax FD-1

Ultrasonic Pocket Fetal Doppler



Compact Battery Charger



Backlit LCD



Built-in Recorder



Interchangeable Waterproof Probe

SonoTrax Series Ultrasonic Pocket Dopplers' superior sensitivity, high-fidelity sound, recording option and optional rechargeable battery ensure the doppler family to meet the routine examination requirements of obstetricians and midwives.

Specifications

- Display: 45 x 25 mm LCD display
- Heart Rate Range: 50~210 bpm
- Resolution: 1 bpm
- Recording Time: 240 s
- Audio Output Power: 1W
- Auto Power-off: 1 minute with no signal or operation
- Size: 34 (D) x 89 (W) x 141 (H) mm
- Weight: about 300 g (including the battery)



Carry Case



Rolling Stand

- Easy-to-use
- Compact design
- Support 1.5 V AA battery
- Built-in audio recorder
- Large backlit LCD display
- Accurate fetal heart detection with high-fidelity sound
- Support 2 MHz, 3 MHz, 4 MHz, 5 MHz, 8 MHz waterproof probes
- Automatic and manual counting of FHR
- Interchangeable waterproof probe
- Automatic power-off



SonoTrax Vascular



SonoTrax Lite



SonoTrax Basic



SonoTrax Basic A



SonoTrax Pro



SonoTrax II



SonoTrax II Pro