

ENDO



Comfort RA660

Front Protection Apron

<https://endo.id>


Ergonomic RA660

Front Protection Apron

Our RA660 range offers up-to-date radiation protection clothing which accompanies your active working day and provides outstanding service. The extraordinarily good wearing comfort is supported by the wide stretch insert integrated in the back panels.

Protection/Lead Equivalent Standard Collection


NovaLite, Leadfree or Standard lead



Front:
0.25 mm Pb, 0.35 mm Pb
or 0.50 mm Pb
Back generally:
0.25 mm Pb

Protection/Lead Equivalent Levio Edition

NovaLite, Leadfree or Standard lead



Total value Levio-Panel:
0.42 mm Pb or 0.50 mm Pb
(outside 0.25 mm Pb)
Back generally:
0.25 mm Pb



- Suitable for a wide range of applications with only frontal radiation exposure
- Ideal fit, maximum comfort for the wearer
- Suitable for long wearing times
- Optionally "without protection" in the back, for a reduced weight of the apron
- Standard equipment includes shoulder pads
- Type testing and certification in accordance with (EU) 2016/425, IEC 61331-1:2014 / DIN EN 61331-1:2016 and IEC 61331-3:2014 / DIN EN 61331-3:2016



Sizes*, recommended for:

- XS Chest girth \leq 84 cm
- S Chest girth approx. 85 - 95 cm
- M Chest girth approx. 96 - 106 cm
- L Chest girth approx. 107 - 117 cm
- XL Chest girth approx. 118 - 128 cm
- XXL Chest girth approx. 129 - 139 cm
- 3XL Chest girth approx. 140 - 150 cm

- Waist/hip girth \leq 91 cm
- Waist/hip girth approx. 92 - 102 cm
- Waist/hip girth approx. 103 - 112 cm
- Waist/hip girth approx. 113 - 123 cm
- Waist/hip girth approx. 124 - 134 cm
- Waist/hip girth approx. 135 - 145 cm
- Waist/hip girth approx. 146 - 156 cm

*The size depends on the girth of the chest, waist, or hip, depending on which of the three measurements is the largest.

Sales Champion

Our customers report that this model is their favourite front apron. This is because it can be put on and taken off easily, fits the body well, stays in place well and provides excellent shoulder relief.



Outlast Temperature Regulation


Closure System

ENDO uses only top quality, low-wear hook-and-loop fasteners.



Stretch insert

The stretch insert running across the entire width of the back ensures a good fit and shoulder relief even when the X-ray protective clothing is worn for long periods.



Higher Protection - Levio Edition

The area with the reinforced protective material is clearly recognisable by its colour.



User Friendly

The sophisticated design of the RA660 front protection apron makes it easy and quick to put on or take off.

Special Feature

The stretch insert offers ideal relief especially for long periods of wear.

Outer Material

Colours
(see front page)



Optional:
ComforTex® HPMF,
ComforTex® HPMF Hybrid

Sizes

From XS to 3XL
(see front page)



Lengths / Body height:
90 cm / 150 - 160 cm
100 cm / 160 - 170 cm
110 cm / 170 - 180 cm
120 cm / 175 - 190 cm
130 cm / 185 - 195 cm
140 cm / 195 - 205 cm

Optional




On request:
Outlast
Breast pocket
Embroidery
Made-to-measure
protection

Accessories

Bolero
Thyroid/sternum



ENDO

 Jalan Raya Menganti 14 Kedurus, Surabaya 60223
Jawa Timur, Indonesia

 62-31-7673636

 62-31-7673737

 0800 . 177 . ENDO (3636)

 info@endo.id

 <https://endo.id>