

ENDO



ENDO

LED Examination Lamp

EI.OL3

<https://endo.id>

ENDO LED Examination Lamp EI.OL3

ENDO LED Examination Lamp EI.OL3 is suitable for exam room.

Specifications

Mode	ENDO EI.OL3
• Illumination at 1 m(lux)	≥80000
• Optional color temperature(K)	6700≥Tc≥3500K
• Input power(VA)	20
• Color rendering index(Ra)	100≥Ra≥85
• Tilt left and right	160°
• Tilt front and back	190°
• Spring arm adjusted	≥600

Features

- Adopt imported LED cold light source
- Adopt digital means for stepless regulation of LED brightness.
- Adjust color temperature and with R9 function
- Unique balance arm design, which is flexible up and down





ENDO EI.OL3




ENDO


 Jalan Raya Menganti 14 Kedurus, Surabaya 60223
Jawa Timur, Indonesia

 62-31-7673636

 info@endo.id

 62-31-7673737

 <https://endo.id>

 0800 . 177 . ENDO (3636)