


ENDO



ENDO

LED Examination Lamp

EI.OL7

<https://endo.id>

ENDO LED Examination Lamp EI.OL7


ENDO LED Examination Lamp EI.OL7 is suitable for exam room.


Specifications

Mode	ENDO EI.OL7
• Illumination at 1 m(lux)	80000
• Color temperature(K)	4000±500
• Color rendering index(Ra)	100≥Ra≥85
• Depth of light beam(cm)	≥500
• The value of total irradiance	<6mW/m2lx
• Total irradiance (Temperature rise on the head of patient under surgical operation ≤ 2°C)	≤600W/m2
• Power supply	220±10%50Hz
• Rated voltage of bulb	24V
• Rated power of bulb	125 W
• Input power(VA)	150





ENDO


 Jalan Raya Menganti 14 Kedurus, Surabaya 60223
Jawa Timur, Indonesia

 62-31-7673636

 info@endo.id

 62-31-7673737

 <https://endo.id>

 0800 . 177 . ENDO (3636)