


# ENDO PO-2

## Finger Oximeter



# ENDO PO-2

## Finger Oximeter


ENDO PO-2 Finger Oximeter delivers a fast, reliable oximetry in a very handy, finger-tip solution by combining the monitor and the sensor into one unit. ENDO PO-2 is small in size, nonetheless it provides SpO<sub>2</sub>, pulse rate measurement, bar graph and plethysmogram display.

- Compact & lightweight design
- Dual color OLED display
- SpO<sub>2</sub>, PR, Bar Graph and Plethysmogram display
- One-button operation
- 4-direction display with 6 display modes
- 10-level adjustable brightness
- Low power consumption
- Automatic power off design
- Low-battery indicator
- Standard 2 AAA batteries power supply



4 Colors

### Features



6 Display Modes



Pocket-size Design



One-button Operation