

ENDO

ENDO Power Wheelchair ENDO JRWD6012



PARAMETERS

Driving range > 20 km

Motor 250W * 2 pcs (Jerry)

Battery 24V 12Ah/20Ah (Lithium Battery)

Parking brake electromagnetic brake

Max speed 6 km/h

Frame material aluminium alloy

Frame color silver, gray

Cushion color black

Control Jerry

Features

Foldable directly

Removable battery box for recharging

Controller can be exchanged from right side to the left

Long driving distance, up to 25 km

Series

ENDO JRWD6012 Standard

ENDO JRWD6012 with Remote Control

ENDO JRWD6012 with Backrest Adjustment