

Generator



Product Code	Compatible with	Packaging
EG GN15	All EG ENERGY® device	1 pcs/case

Shear Scalpel



Product Code	Shaft Diameter	Shaft Length	Tip	Packaging
EG US 36D	5mm	36cm	Curved	3 pcs/case
EG US 23D	5mm	23cm	Curved	3 pcs/case
EG US 14D	5mm	14cm	Curved	3 pcs/case



Product Code	Shaft Diameter	Shaft Length	Tip	Packaging
EG US 36R	5mm	36cm	Curved	3 pcs/case
EG US 23R	5mm	23cm	Curved	3 pcs/case
EG US 14R	5mm	14cm	Curved	3 pcs/case

Forceps Scalpel



Product Code	Length	Tip	Packaging
EG UF 09D	9mm	Curved	5 pcs/case



Product Code	Length	Tip	Packaging
EG UF 09R	9mm	Curved	5 pcs/case

Transducer



Product Code	Compatible with	Packaging
EG TR100	EG US 36D, EG US 23D, EG US 14D EG US 36R, EG US 23R, EG US 14R	1 pcs/case



Product Code	Compatible with	Packaging
EG TR200	EG UF 09D, EG UF 09R	1 pcs/case



Engine Medical Equipment Manufacturing (Shanghai) Corporation.

Add : No.1199, Xinsiping Road, 201413, Shanghai, China


Tel : +86-21-57125335-8008 Fax : +86-21-57125339

Https://www.enginemed.com



EG ENERGY™
Ultrasonic Surgical System

EG ENERGY™ Ultrasonic Surgical System



Protection Pad
Effectively extend the using time of the scalpel

Integrated Knob
Equipped with torque wrench function. Easy to connect the scalpel to transducer

Precise Dissection
Lateral thermal tissue damage: Only 1mm

Multifuction
Grasping, dissection, cutting, coagulation.....All in one single instrument

First Forceps made in china


Faster cutting speed

Better coagulation

More compatible




Faster cutting speed



Average cutting speed : 4.3s(EG US36) Vs. 5.4s(similar product)

Reliably sealing vessels up to 5 mm



Burst pressure for 3-5mm vessels: EG US 36 vs. similar product.