




## Oscillating Saw Blades for Orthopaedics




## Reciprocating Saw Blades for Bone Surgery



## Saw Blades for Trauma



## Saw Blades for small Bones



MANUFACTURER

DISTRIBUTOR

**GOMINA**  
SWISS QUALITY

**GOMINA AG**  
Raiftstrasse 4  
CH-3989 Niederwald  
Telephone ++41 (0)27 970 12 70  
Telefax ++41 (0)27 970 12 75  
[info@gomina.ch](mailto:info@gomina.ch)  
[www.gomina.ch](http://www.gomina.ch)

11/2017

## Table of Content

Gomina Saw Blades .....	3
Premium Line .....	4
Saw Blade fits to Aesculap®.....	5
Saw Blade fits to Linvatec Hall®.....	9
Saw Blade fits to AO (Synthes®, Sodem®).....	13
Saw Blade fits to deSoutter (Multidrive™).....	17
Saw Blade fits to Stryker®, System EHD, 2000 / 4 / 5 / 6 / 7.....	19
Saw Blade fits to Aesculap® Acculan 3Ti.....	24
Saw Blade fits to MicroAire® (Series 7000).....	28
Saw Blade fits to Arthrex Power System 600 .....	30
Performance Line .....	31
Saw Blade fits to Aesculap®.....	32
Saw Blade fits to Linvatec Hall®.....	34
Saw Blade fits to AO (Synthes®, Sodem®).....	37
Saw Blade fits to deSoutter (Multidrive™).....	41
Saw Blade fits to Stryker®, System EHD, 2000 / 4 / 5 / 6.....	44
Saw Blade fits to Aesculap® Acculan 3Ti.....	48
Saw Blade fits to MicroAire® (Series 7000).....	51
Reciprocating Saw Blades.....	53
Reciprocating Saw Blade fits to Linvatec Hall®, Power Pro, Versipower .....	54
Reciprocating Saw Blade fits to Stryker®, System EHD, 2000 / 4 / 5 / 6 / 7.....	55
Reciprocating Saw Blade fits to Tuke® .....	56
Reciprocating Saw Blade fits to Synthes®.....	57
Reciprocating Saw Blade fits to Zimmer® .....	60
Reciprocating Saw Blade fits to Aesculap® Acculan 3Ti .....	61
Reciprocating Keel Blade .....	62
Reciprocating Keel Blade fits to Linvatec Hall®, Power Pro, Versipower .....	63
Reciprocating Keel Blade fits to Stryker®, System EHD, 2000 / 4 / 5 / 6 / 7.....	64
Reciprocating Keel Blade fits to Synthes® .....	65
Reciprocating Keel Blade fits to Aesculap® Acculan 3Ti.....	66
Reciprocating Keel Blade fits to Zimmer® .....	67
Saw Blades for Trauma .....	68
Saw Blade fits to Aesculap®.....	69
Saw Blade fits to Linvatec Hall®.....	71
Saw Blade fits to AO (Synthes®, Sodem®).....	72
Saw Blade fits to deSoutter (Multidrive™).....	76
Saw Blade fits to Stryker®, System EHD, 2000 / 4 / 5 / 6 / 7.....	77
Saw Blade fits to Aesculap® Acculan 3Ti.....	80
Saw Blades for small Bones .....	83
Saw Blade fits to Stryker®, TPS.....	84
Saw Blade fits to Smith & Nephew DYONICS™ Saw REF 7205786 and 7205791 .....	86
Saw Blade fits to Linvatec Hall®, MicroChoice®, Micro100™, Micro E® .....	87
Saw Blade fits to Colibri (Synthes®).....	89
Saw Blade fits to EPEN (Synthes®) .....	92
Saw Blade fits to Arthrex .....	94
Saw Blade fits to Aesculap® Acculan 3Ti.....	95

## Gomina Saw Blades

### NON STERILE



### STERILE

Saw blade, double sterile packaged, with teeth protection, in boxes of 5 pieces.  
6 folded patient labels for each blade, placed between first and second pouch.



## Premium Line

Manufacturer: **GOMINA**

Cutting Thickness

Cutting length in mm




Swiss Quality

Reference-No


LOT-Nr.

Gomina Teeth-  
Design Premium Line


## HUB OVERVIEW




Aesculap®




Stryker®, System EHD,  
2000 / 4 / 5 / 6 / 7




Linvatec Hall®




Aesculap®  
Acculan 3Ti




AO (Synthes®,  
Soderm®)



MicroAire®  
(Series 7000)



deSoutter  
(Multidrive™)



Arthrex Power  
System 600

## Saw Blade fits to Aesculap®

**90 MM LENGTH      25 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
90 mm	25 mm	<b>0.90 mm</b>	<b>212.259.090</b>	<b>212.259.090S</b>
90 mm	25 mm	<b>1.00 mm</b>	<b>212.259.100</b>	<b>212.259.100S</b>
90 mm	25 mm	<b>1.19 mm</b>	<b>212.259.119</b>	<b>212.259.119S</b>
90 mm	25 mm	<b>1.23 mm</b>	<b>212.259.123</b>	<b>212.259.123S</b>
90 mm	25 mm	<b>1.25 mm</b>	<b>212.259.125</b>	<b>212.259.125S</b>
90 mm	25 mm	<b>1.27 mm</b>	<b>212.259.127</b>	<b>212.259.127S</b>
90 mm	25 mm	<b>1.37 mm</b>	<b>212.259.137</b>	<b>212.259.137S</b>
90 mm	25 mm	<b>1.47 mm</b>	<b>212.259.147</b>	<b>212.259.147S</b>

**90 MM LENGTH      19 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
90 mm	19 mm	<b>0.90 mm</b>	<b>212.199.090</b>	<b>212.199.090S</b>
90 mm	19 mm	<b>1.00 mm</b>	<b>212.199.100</b>	<b>212.199.100S</b>
90 mm	19 mm	<b>1.19 mm</b>	<b>212.199.119</b>	<b>212.199.119S</b>
90 mm	19 mm	<b>1.23 mm</b>	<b>212.199.123</b>	<b>212.199.123S</b>
90 mm	19 mm	<b>1.25 mm</b>	<b>212.199.125</b>	<b>212.199.125S</b>
90 mm	19 mm	<b>1.27 mm</b>	<b>212.199.127</b>	<b>212.199.127S</b>
90 mm	19 mm	<b>1.37 mm</b>	<b>212.199.137</b>	<b>212.199.137S</b>
90 mm	19 mm	<b>1.47 mm</b>	<b>212.199.147</b>	<b>212.199.147S</b>

## Saw Blade fits to Aesculap®

**90 MM LENGTH      12 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
90 mm	12 mm	<b>0.90 mm</b>	<b>212.129.090</b>	<b>212.129.090S</b>
90 mm	12 mm	<b>1.00 mm</b>	<b>212.129.100</b>	<b>212.129.100S</b>
90 mm	12 mm	<b>1.19 mm</b>	<b>212.129.119</b>	<b>212.129.119S</b>
90 mm	12 mm	<b>1.23 mm</b>	<b>212.129.123</b>	<b>212.129.123S</b>
90 mm	12 mm	<b>1.25 mm</b>	<b>212.129.125</b>	<b>212.129.125S</b>
90 mm	12 mm	<b>1.27 mm</b>	<b>212.129.127</b>	<b>212.129.127S</b>
90 mm	12 mm	<b>1.37 mm</b>	<b>212.129.137</b>	<b>212.129.137S</b>
90 mm	12 mm	<b>1.47 mm</b>	<b>212.129.147</b>	<b>212.129.147S</b>

**90 MM LENGTH      12 - 19 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
90 mm	12 - 19 mm	<b>0.90 mm</b>	<b>212.999.090</b>	<b>212.999.090S</b>
90 mm	12 - 19 mm	<b>1.00 mm</b>	<b>212.999.100</b>	<b>212.999.100S</b>
90 mm	12 - 19 mm	<b>1.19 mm</b>	<b>212.999.119</b>	<b>212.999.119S</b>
90 mm	12 - 19 mm	<b>1.23 mm</b>	<b>212.999.123</b>	<b>212.999.123S</b>
90 mm	12 - 19 mm	<b>1.25 mm</b>	<b>212.999.125</b>	<b>212.999.125S</b>
90 mm	12 - 19 mm	<b>1.27 mm</b>	<b>212.999.127</b>	<b>212.999.127S</b>
90 mm	12 - 19 mm	<b>1.37 mm</b>	<b>212.999.137</b>	<b>212.999.137S</b>
90 mm	12 - 19 mm	<b>1.47 mm</b>	<b>212.999.147</b>	<b>212.999.147S</b>

## Saw Blade fits to Aesculap®

**70 MM LENGTH      25 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
70 mm	25 mm	<b>0.90 mm</b>	<b>212.257.090</b>	<b>212.257.090S</b>
70 mm	25 mm	<b>1.00 mm</b>	<b>212.257.100</b>	<b>212.257.100S</b>
70 mm	25 mm	<b>1.19 mm</b>	<b>212.257.119</b>	<b>212.257.119S</b>
70 mm	25 mm	<b>1.23 mm</b>	<b>212.257.123</b>	<b>212.257.123S</b>
70 mm	25 mm	<b>1.25 mm</b>	<b>212.257.125</b>	<b>212.257.125S</b>
70 mm	25 mm	<b>1.27 mm</b>	<b>212.257.127</b>	<b>212.257.127S</b>
70 mm	25 mm	<b>1.37 mm</b>	<b>212.257.137</b>	<b>212.257.137S</b>
70 mm	25 mm	<b>1.47 mm</b>	<b>212.257.147</b>	<b>212.257.147S</b>

**70 MM LENGTH      19 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
70 mm	19 mm	<b>0.90 mm</b>	<b>212.197.090</b>	<b>212.197.090S</b>
70 mm	19 mm	<b>1.00 mm</b>	<b>212.197.100</b>	<b>212.197.100S</b>
70 mm	19 mm	<b>1.19 mm</b>	<b>212.197.119</b>	<b>212.197.119S</b>
70 mm	19 mm	<b>1.23 mm</b>	<b>212.197.123</b>	<b>212.197.123S</b>
70 mm	19 mm	<b>1.25 mm</b>	<b>212.197.125</b>	<b>212.197.125S</b>
70 mm	19 mm	<b>1.27 mm</b>	<b>212.197.127</b>	<b>212.197.127S</b>
70 mm	19 mm	<b>1.37 mm</b>	<b>212.197.137</b>	<b>212.197.137S</b>
70 mm	19 mm	<b>1.47 mm</b>	<b>212.197.147</b>	<b>212.197.147S</b>

## Saw Blade fits to Aesculap®

50 MM LENGTH

25 MM WIDTH



Length	Width	Thickness	REF No non-sterile	REF No sterile
50 mm	25 mm	<b>0.90 mm</b>	<b>212.255.090</b>	<b>212.255.090S</b>



## Saw Blade fits to Linvatec Hall®

**100 MM LENGTH    19 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
100 mm	19 mm	1.27 mm	213.190.127	213.190.127S

**90 MM LENGTH    25 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
90 mm	25 mm	0.90 mm	213.259.090	213.259.090S
90 mm	25 mm	1.00 mm	213.259.100	213.259.100S
90 mm	25 mm	1.19 mm	213.259.119	213.259.119S
90 mm	25 mm	1.23 mm	213.259.123	213.259.123S
90 mm	25 mm	1.25 mm	213.259.125	213.259.125S
90 mm	25 mm	1.27 mm	213.259.127	213.259.127S
90 mm	25 mm	1.37 mm	213.259.137	213.259.137S
90 mm	25 mm	1.47 mm	213.259.147	213.259.147S

## Saw Blade fits to Linvatec Hall®

**90 MM LENGTH      19 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
90 mm	19 mm	<b>0.90 mm</b>	<b>213.199.090</b>	<b>213.199.090S</b>
90 mm	19 mm	<b>1.00 mm</b>	<b>213.199.100</b>	<b>213.199.100S</b>
90 mm	19 mm	<b>1.19 mm</b>	<b>213.199.119</b>	<b>213.199.119S</b>
90 mm	19 mm	<b>1.23 mm</b>	<b>213.199.123</b>	<b>213.199.123S</b>
90 mm	19 mm	<b>1.25 mm</b>	<b>213.199.125</b>	<b>213.199.125S</b>
90 mm	19 mm	<b>1.27 mm</b>	<b>213.199.127</b>	<b>213.199.127S</b>
90 mm	19 mm	<b>1.37 mm</b>	<b>213.199.137</b>	<b>213.199.137S</b>
90 mm	19 mm	<b>1.47 mm</b>	<b>213.199.147</b>	<b>213.199.147S</b>

**90 MM LENGTH      12 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
90 mm	12 mm	<b>0.90 mm</b>	<b>213.129.090</b>	<b>213.129.090S</b>
90 mm	12 mm	<b>1.00 mm</b>	<b>213.129.100</b>	<b>213.129.100S</b>
90 mm	12 mm	<b>1.19 mm</b>	<b>213.129.119</b>	<b>213.129.119S</b>
90 mm	12 mm	<b>1.23 mm</b>	<b>213.129.123</b>	<b>213.129.123S</b>
90 mm	12 mm	<b>1.25 mm</b>	<b>213.129.125</b>	<b>213.129.125S</b>
90 mm	12 mm	<b>1.27 mm</b>	<b>213.129.127</b>	<b>213.129.127S</b>
90 mm	12 mm	<b>1.37 mm</b>	<b>213.129.137</b>	<b>213.129.137S</b>
90 mm	12 mm	<b>1.47 mm</b>	<b>213.129.147</b>	<b>213.129.147S</b>

## Saw Blade fits to Linvatec Hall®

**90 MM LENGTH      12 - 19 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
90 mm	12 - 19 mm	<b>0.90 mm</b>	<b>213.999.090</b>	<b>213.999.090S</b>
90 mm	12 - 19 mm	<b>1.00 mm</b>	<b>213.999.100</b>	<b>213.999.100S</b>
90 mm	12 - 19 mm	<b>1.19 mm</b>	<b>213.999.119</b>	<b>213.999.119S</b>
90 mm	12 - 19 mm	<b>1.23 mm</b>	<b>213.999.123</b>	<b>213.999.123S</b>
90 mm	12 - 19 mm	<b>1.25 mm</b>	<b>213.999.125</b>	<b>213.999.125S</b>
90 mm	12 - 19 mm	<b>1.27 mm</b>	<b>213.999.127</b>	<b>213.999.127S</b>
90 mm	12 - 19 mm	<b>1.37 mm</b>	<b>213.999.137</b>	<b>213.999.137S</b>
90 mm	12 - 19 mm	<b>1.47 mm</b>	<b>213.999.147</b>	<b>213.999.147S</b>

**70 MM LENGTH      25 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
70 mm	25 mm	<b>0.90 mm</b>	<b>213.257.090</b>	<b>213.257.090S</b>
70 mm	25 mm	<b>1.00 mm</b>	<b>213.257.100</b>	<b>213.257.100S</b>
70 mm	25 mm	<b>1.19 mm</b>	<b>213.257.119</b>	<b>213.257.119S</b>
70 mm	25 mm	<b>1.23 mm</b>	<b>213.257.123</b>	<b>213.257.123S</b>
70 mm	25 mm	<b>1.25 mm</b>	<b>213.257.125</b>	<b>213.257.125S</b>
70 mm	25 mm	<b>1.27 mm</b>	<b>213.257.127</b>	<b>213.257.127S</b>
70 mm	25 mm	<b>1.37 mm</b>	<b>213.257.137</b>	<b>213.257.137S</b>
70 mm	25 mm	<b>1.47 mm</b>	<b>213.257.147</b>	<b>213.257.147S</b>

## Saw Blade fits to Linvatec Hall®

**70 MM LENGTH      19 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
70 mm	19 mm	<b>0.90 mm</b>	<b>213.197.090</b>	<b>213.197.090S</b>
70 mm	19 mm	<b>1.00 mm</b>	<b>213.197.100</b>	<b>213.197.100S</b>
70 mm	19 mm	<b>1.19 mm</b>	<b>213.197.119</b>	<b>213.197.119S</b>
70 mm	19 mm	<b>1.23 mm</b>	<b>213.197.123</b>	<b>213.197.123S</b>
70 mm	19 mm	<b>1.25 mm</b>	<b>213.197.125</b>	<b>213.197.125S</b>
70 mm	19 mm	<b>1.27 mm</b>	<b>213.197.127</b>	<b>213.197.127S</b>
70 mm	19 mm	<b>1.37 mm</b>	<b>213.197.137</b>	<b>213.197.137S</b>
70 mm	19 mm	<b>1.47 mm</b>	<b>213.197.147</b>	<b>213.197.147S</b>

**70 MM LENGTH      12 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
70 mm	12 mm	<b>1.00 mm</b>	<b>213.127.100</b>	<b>213.127.100S</b>


## Saw Blade fits to AO (Synthes®, Sodem®)

**90 MM LENGTH      25 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
90 mm	25 mm	<b>0.90 mm</b>	<b>215.259.090</b>	<b>215.259.090S</b>
90 mm	25 mm	<b>1.00 mm</b>	<b>215.259.100</b>	<b>215.259.100S</b>
90 mm	25 mm	<b>1.19 mm</b>	<b>215.259.119</b>	<b>215.259.119S</b>
90 mm	25 mm	<b>1.23 mm</b>	<b>215.259.123</b>	<b>215.259.123S</b>
90 mm	25 mm	<b>1.25 mm</b>	<b>215.259.125</b>	<b>215.259.125S</b>
90 mm	25 mm	<b>1.27 mm</b>	<b>215.259.127</b>	<b>215.259.127S</b>
90 mm	25 mm	<b>1.37 mm</b>	<b>215.259.137</b>	<b>215.259.137S</b>
90 mm	25 mm	<b>1.47 mm</b>	<b>215.259.147</b>	<b>215.259.147S</b>

**90 MM LENGTH      22 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
90 mm	22 mm	<b>1.19 mm</b>	<b>215.229.119</b>	<b>215.229.119S</b>
90 mm	22 mm	<b>1.27 mm</b>	<b>215.229.127</b>	<b>215.229.127S</b>

## Saw Blade fits to AO (Synthes®, Sodem®)

**90 MM LENGTH      19 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
90 mm	19 mm	<b>0.90 mm</b>	<b>215.199.090</b>	<b>215.199.090S</b>
90 mm	19 mm	<b>1.00 mm</b>	<b>215.199.100</b>	<b>215.199.100S</b>
90 mm	19 mm	<b>1.19 mm</b>	<b>215.199.119</b>	<b>215.199.119S</b>
90 mm	19 mm	<b>1.23 mm</b>	<b>215.199.123</b>	<b>215.199.123S</b>
90 mm	19 mm	<b>1.25 mm</b>	<b>215.199.125</b>	<b>215.199.125S</b>
90 mm	19 mm	<b>1.27 mm</b>	<b>215.199.127</b>	<b>215.199.127S</b>
90 mm	19 mm	<b>1.37 mm</b>	<b>215.199.137</b>	<b>215.199.137S</b>
90 mm	19 mm	<b>1.47 mm</b>	<b>215.199.147</b>	<b>215.199.147S</b>

**90 MM LENGTH      12 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
90 mm	12 mm	<b>0.90 mm</b>	<b>215.129.090</b>	<b>215.129.090S</b>
90 mm	12 mm	<b>1.00 mm</b>	<b>215.129.100</b>	<b>215.129.100S</b>
90 mm	12 mm	<b>1.19 mm</b>	<b>215.129.119</b>	<b>215.129.119S</b>
90 mm	12 mm	<b>1.23 mm</b>	<b>215.129.123</b>	<b>215.129.123S</b>
90 mm	12 mm	<b>1.25 mm</b>	<b>215.129.125</b>	<b>215.129.125S</b>
90 mm	12 mm	<b>1.27 mm</b>	<b>215.129.127</b>	<b>215.129.127S</b>
90 mm	12 mm	<b>1.37 mm</b>	<b>215.129.137</b>	<b>215.129.137S</b>
90 mm	12 mm	<b>1.47 mm</b>	<b>215.129.147</b>	<b>215.129.147S</b>

## Saw Blade fits to AO (Synthes®, Sodem®)

**90 MM LENGTH      12 - 19 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
90 mm	12 - 19 mm	<b>0.90 mm</b>	<b>215.999.090</b>	<b>215.999.090S</b>
90 mm	12 - 19 mm	<b>1.00 mm</b>	<b>215.999.100</b>	<b>215.999.100S</b>
90 mm	12 - 19 mm	<b>1.19 mm</b>	<b>215.999.119</b>	<b>215.999.119S</b>
90 mm	12 - 19 mm	<b>1.23 mm</b>	<b>215.999.123</b>	<b>215.999.123S</b>
90 mm	12 - 19 mm	<b>1.25 mm</b>	<b>215.999.125</b>	<b>215.999.125S</b>
90 mm	12 - 19 mm	<b>1.27 mm</b>	<b>215.999.127</b>	<b>215.999.127S</b>
90 mm	12 - 19 mm	<b>1.37 mm</b>	<b>215.999.137</b>	<b>215.999.137S</b>
90 mm	12 - 19 mm	<b>1.47 mm</b>	<b>215.999.147</b>	<b>215.999.147S</b>

**70 MM LENGTH      25 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
70 mm	25 mm	<b>0.90 mm</b>	<b>215.257.090</b>	<b>215.257.090S</b>
70 mm	25 mm	<b>1.00 mm</b>	<b>215.257.100</b>	<b>215.257.100S</b>
70 mm	25 mm	<b>1.19 mm</b>	<b>215.257.119</b>	<b>215.257.119S</b>
70 mm	25 mm	<b>1.23 mm</b>	<b>215.257.123</b>	<b>215.257.123S</b>
70 mm	25 mm	<b>1.25 mm</b>	<b>215.257.125</b>	<b>215.257.125S</b>
70 mm	25 mm	<b>1.27 mm</b>	<b>215.257.127</b>	<b>215.257.127S</b>
70 mm	25 mm	<b>1.37 mm</b>	<b>215.257.137</b>	<b>215.257.137S</b>
70 mm	25 mm	<b>1.47 mm</b>	<b>215.257.147</b>	<b>215.257.147S</b>

## Saw Blade fits to AO (Synthes®, Sodem®)

**70 MM LENGTH      19 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
70 mm	19 mm	0.90 mm	215.197.090	215.197.090S
70 mm	19 mm	1.00 mm	215.197.100	215.197.100S
70 mm	19 mm	1.19 mm	215.197.119	215.197.119S
70 mm	19 mm	1.23 mm	215.197.123	215.197.123S
70 mm	19 mm	1.25 mm	215.197.125	215.197.125S
70 mm	19 mm	1.27 mm	215.197.127	215.197.127S
70 mm	19 mm	1.37 mm	215.197.137	215.197.137S
70 mm	19 mm	1.47 mm	215.197.147	215.197.147S

**70 MM LENGTH      6 - 12 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
70 mm	6 - 12 mm	1.27 mm	215.987.127	215.987.127S

**50 MM LENGTH      25 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
50 mm	25 mm	0.90 mm	215.255.090	215.255.090S



## Saw Blade fits to deSoutter (Multidrive™)

**90 MM LENGTH      25 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
90 mm	25 mm	<b>0.90 mm</b>	<b>217.259.090</b>	<b>217.259.090S</b>
90 mm	25 mm	<b>1.00 mm</b>	<b>217.259.100</b>	<b>217.259.100S</b>
90 mm	25 mm	<b>1.19 mm</b>	<b>217.259.119</b>	<b>217.259.119S</b>
90 mm	25 mm	<b>1.23 mm</b>	<b>217.259.123</b>	<b>217.259.123S</b>
90 mm	25 mm	<b>1.25 mm</b>	<b>217.259.125</b>	<b>217.259.125S</b>
90 mm	25 mm	<b>1.27 mm</b>	<b>217.259.127</b>	<b>217.259.127S</b>
90 mm	25 mm	<b>1.37 mm</b>	<b>217.259.137</b>	<b>217.259.137S</b>
90 mm	25 mm	<b>1.47 mm</b>	<b>217.259.147</b>	<b>217.259.147S</b>

**90 MM LENGTH      19 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
90 mm	19 mm	<b>0.90 mm</b>	<b>217.199.090</b>	<b>217.199.090S</b>
90 mm	19 mm	<b>1.00 mm</b>	<b>217.199.100</b>	<b>217.199.100S</b>
90 mm	19 mm	<b>1.19 mm</b>	<b>217.199.119</b>	<b>217.199.119S</b>
90 mm	19 mm	<b>1.23 mm</b>	<b>217.199.123</b>	<b>217.199.123S</b>
90 mm	19 mm	<b>1.25 mm</b>	<b>217.199.125</b>	<b>217.199.125S</b>
90 mm	19 mm	<b>1.27 mm</b>	<b>217.199.127</b>	<b>217.199.127S</b>
90 mm	19 mm	<b>1.37 mm</b>	<b>217.199.137</b>	<b>217.199.137S</b>
90 mm	19 mm	<b>1.47 mm</b>	<b>217.199.147</b>	<b>217.199.147S</b>

## Saw Blade fits to deSoutter (Multidrive™)

**90 MM LENGTH      12 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
90 mm	12 mm	<b>0.90 mm</b>	<b>217.129.090</b>	<b>217.129.090S</b>
90 mm	12 mm	<b>1.00 mm</b>	<b>217.129.100</b>	<b>217.129.100S</b>
90 mm	12 mm	<b>1.19 mm</b>	<b>217.129.119</b>	<b>217.129.119S</b>
90 mm	12 mm	<b>1.23 mm</b>	<b>217.129.123</b>	<b>217.129.123S</b>
90 mm	12 mm	<b>1.25 mm</b>	<b>217.129.125</b>	<b>217.129.125S</b>
90 mm	12 mm	<b>1.27 mm</b>	<b>217.129.127</b>	<b>217.129.127S</b>
90 mm	12 mm	<b>1.37 mm</b>	<b>217.129.137</b>	<b>217.129.137S</b>
90 mm	12 mm	<b>1.47 mm</b>	<b>217.129.147</b>	<b>217.129.147S</b>


**90 MM LENGTH      12 - 19 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
90 mm	12 - 19 mm	<b>0.90 mm</b>	<b>217.999.090</b>	<b>217.999.090S</b>
90 mm	12 - 19 mm	<b>1.00 mm</b>	<b>217.999.100</b>	<b>217.999.100S</b>
90 mm	12 - 19 mm	<b>1.19 mm</b>	<b>217.999.119</b>	<b>217.999.119S</b>
90 mm	12 - 19 mm	<b>1.23 mm</b>	<b>217.999.123</b>	<b>217.999.123S</b>
90 mm	12 - 19 mm	<b>1.25 mm</b>	<b>217.999.125</b>	<b>217.999.125S</b>
90 mm	12 - 19 mm	<b>1.27 mm</b>	<b>217.999.127</b>	<b>217.999.127S</b>
90 mm	12 - 19 mm	<b>1.37 mm</b>	<b>217.999.137</b>	<b>217.999.137S</b>
90 mm	12 - 19 mm	<b>1.47 mm</b>	<b>217.999.147</b>	<b>217.999.147S</b>


**Saw Blade fits to Stryker®, System EHD, 2000 / 4 / 5 / 6 / 7**

**110 MM LENGTH    25 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
110 mm	25 mm	1.27 mm	218.251.127	218.251.127S

**110 MM LENGTH    19 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
110 mm	19 mm	1.27 mm	218.191.127	218.191.127S

**100 MM LENGTH    25 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
100 mm	25 mm	1.19 mm	218.250.119	218.250.119S


**100 MM LENGTH    19 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
100 mm	19 mm	1.27 mm	218.190.127	218.190.127S


## Saw Blade fits to Stryker®, System EHD, 2000 / 4 / 5 / 6 / 7

**90 MM LENGTH      25 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
90 mm	25 mm	<b>0.90 mm</b>	<b>218.259.090</b>	<b>218.259.090S</b>
90 mm	25 mm	<b>1.00 mm</b>	<b>218.259.100</b>	<b>218.259.100S</b>
90 mm	25 mm	<b>1.19 mm</b>	<b>218.259.119</b>	<b>218.259.119S</b>
90 mm	25 mm	<b>1.23 mm</b>	<b>218.259.123</b>	<b>218.259.123S</b>
90 mm	25 mm	<b>1.25 mm</b>	<b>218.259.125</b>	<b>218.259.125S</b>
90 mm	25 mm	<b>1.27 mm</b>	<b>218.259.127</b>	<b>218.259.127S</b>
90 mm	25 mm	<b>1.37 mm</b>	<b>218.259.137</b>	<b>218.259.137S</b>
90 mm	25 mm	<b>1.47 mm</b>	<b>218.259.147</b>	<b>218.259.147S</b>


**90 MM LENGTH      22 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
90 mm	22 mm	<b>1.27 mm</b>	<b>218.229.127</b>	<b>218.229.127S</b>


## Saw Blade fits to Stryker®, System EHD, 2000 / 4 / 5 / 6 / 7

**90 MM LENGTH      19 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
90 mm	19 mm	<b>0.90 mm</b>	<b>218.199.090</b>	<b>218.199.090S</b>
90 mm	19 mm	<b>1.00 mm</b>	<b>218.199.100</b>	<b>218.199.100S</b>
90 mm	19 mm	<b>1.19 mm</b>	<b>218.199.119</b>	<b>218.199.119S</b>
90 mm	19 mm	<b>1.23 mm</b>	<b>218.199.123</b>	<b>218.199.123S</b>
90 mm	19 mm	<b>1.25 mm</b>	<b>218.199.125</b>	<b>218.199.125S</b>
90 mm	19 mm	<b>1.27 mm</b>	<b>218.199.127</b>	<b>218.199.127S</b>
90 mm	19 mm	<b>1.37 mm</b>	<b>218.199.137</b>	<b>218.199.137S</b>
90 mm	19 mm	<b>1.47 mm</b>	<b>218.199.147</b>	<b>218.199.147S</b>


**90 MM LENGTH      12 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
90 mm	12 mm	<b>0.90 mm</b>	<b>218.129.090</b>	<b>218.129.090S</b>
90 mm	12 mm	<b>1.00 mm</b>	<b>218.129.100</b>	<b>218.129.100S</b>
90 mm	12 mm	<b>1.19 mm</b>	<b>218.129.119</b>	<b>218.129.119S</b>
90 mm	12 mm	<b>1.23 mm</b>	<b>218.129.123</b>	<b>218.129.123S</b>
90 mm	12 mm	<b>1.25 mm</b>	<b>218.129.125</b>	<b>218.129.125S</b>
90 mm	12 mm	<b>1.27 mm</b>	<b>218.129.127</b>	<b>218.129.127S</b>
90 mm	12 mm	<b>1.37 mm</b>	<b>218.129.137</b>	<b>218.129.137S</b>
90 mm	12 mm	<b>1.47 mm</b>	<b>218.129.147</b>	<b>218.129.147S</b>

## Saw Blade fits to Stryker®, System EHD, 2000 / 4 / 5 / 6 / 7

**90 MM LENGTH      12 - 19 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
90 mm	12 - 19 mm	<b>0.90 mm</b>	<b>218.999.090</b>	<b>218.999.090S</b>
90 mm	12 - 19 mm	<b>1.00 mm</b>	<b>218.999.100</b>	<b>218.999.100S</b>
90 mm	12 - 19 mm	<b>1.19 mm</b>	<b>218.999.119</b>	<b>218.999.119S</b>
90 mm	12 - 19 mm	<b>1.23 mm</b>	<b>218.999.123</b>	<b>218.999.123S</b>
90 mm	12 - 19 mm	<b>1.25 mm</b>	<b>218.999.125</b>	<b>218.999.125S</b>
90 mm	12 - 19 mm	<b>1.27 mm</b>	<b>218.999.127</b>	<b>218.999.127S</b>
90 mm	12 - 19 mm	<b>1.37 mm</b>	<b>218.999.137</b>	<b>218.999.137S</b>
90 mm	12 - 19 mm	<b>1.47 mm</b>	<b>218.999.147</b>	<b>218.999.147S</b>

**70 MM LENGTH      25 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
70 mm	25 mm	<b>0.90 mm</b>	<b>218.257.090</b>	<b>218.257.090S</b>

## Saw Blade fits to Stryker®, System EHD, 2000 / 4 / 5 / 6 / 7

**70 MM LENGTH      19 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
70 mm	19 mm	<b>0.90 mm</b>	<b>218.197.090</b>	<b>218.197.090S</b>
70 mm	19 mm	<b>1.00 mm</b>	<b>218.197.100</b>	<b>218.197.100S</b>
70 mm	19 mm	<b>1.19 mm</b>	<b>218.197.119</b>	<b>218.197.119S</b>
70 mm	19 mm	<b>1.23 mm</b>	<b>218.197.123</b>	<b>218.197.123S</b>
70 mm	19 mm	<b>1.25 mm</b>	<b>218.197.125</b>	<b>218.197.125S</b>
70 mm	19 mm	<b>1.27 mm</b>	<b>218.197.127</b>	<b>218.197.127S</b>
70 mm	19 mm	<b>1.37 mm</b>	<b>218.197.137</b>	<b>218.197.137S</b>
70 mm	19 mm	<b>1.47 mm</b>	<b>218.197.147</b>	<b>218.197.147S</b>

## Saw Blade fits to Aesculap® Acculan 3Ti

**90 MM LENGTH      25 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
90 mm	25 mm	<b>0.90 mm</b>	<b>219.259.090</b>	<b>219.259.090S</b>
90 mm	25 mm	<b>1.00 mm</b>	<b>219.259.100</b>	<b>219.259.100S</b>
90 mm	25 mm	<b>1.19 mm</b>	<b>219.259.119</b>	<b>219.259.119S</b>
90 mm	25 mm	<b>1.23 mm</b>	<b>219.259.123</b>	<b>219.259.123S</b>
90 mm	25 mm	<b>1.25 mm</b>	<b>219.259.125</b>	<b>219.259.125S</b>
90 mm	25 mm	<b>1.27 mm</b>	<b>219.259.127</b>	<b>219.259.127S</b>
90 mm	25 mm	<b>1.37 mm</b>	<b>219.259.137</b>	<b>219.259.137S</b>
90 mm	25 mm	<b>1.47 mm</b>	<b>219.259.147</b>	<b>219.259.147S</b>

**90 MM LENGTH      19 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
90 mm	19 mm	<b>0.90 mm</b>	<b>219.199.090</b>	<b>219.199.090S</b>
90 mm	19 mm	<b>1.00 mm</b>	<b>219.199.100</b>	<b>219.199.100S</b>
90 mm	19 mm	<b>1.19 mm</b>	<b>219.199.119</b>	<b>219.199.119S</b>
90 mm	19 mm	<b>1.23 mm</b>	<b>219.199.123</b>	<b>219.199.123S</b>
90 mm	19 mm	<b>1.25 mm</b>	<b>219.199.125</b>	<b>219.199.125S</b>
90 mm	19 mm	<b>1.27 mm</b>	<b>219.199.127</b>	<b>219.199.127S</b>
90 mm	19 mm	<b>1.37 mm</b>	<b>219.199.137</b>	<b>219.199.137S</b>
90 mm	19 mm	<b>1.47 mm</b>	<b>219.199.147</b>	<b>219.199.147S</b>



## Saw Blade fits to Aesculap® Acculan 3Ti

**90 MM LENGTH      12 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
90 mm	12 mm	<b>0.90 mm</b>	<b>219.129.090</b>	<b>219.129.090S</b>
90 mm	12 mm	<b>1.00 mm</b>	<b>219.129.100</b>	<b>219.129.100S</b>
90 mm	12 mm	<b>1.19 mm</b>	<b>219.129.119</b>	<b>219.129.119S</b>
90 mm	12 mm	<b>1.23 mm</b>	<b>219.129.123</b>	<b>219.129.123S</b>
90 mm	12 mm	<b>1.25 mm</b>	<b>219.129.125</b>	<b>219.129.125S</b>
90 mm	12 mm	<b>1.27 mm</b>	<b>219.129.127</b>	<b>219.129.127S</b>
90 mm	12 mm	<b>1.37 mm</b>	<b>219.129.137</b>	<b>219.129.137S</b>
90 mm	12 mm	<b>1.47 mm</b>	<b>219.129.147</b>	<b>219.129.147S</b>

**90 MM LENGTH      19 - 22 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
90 mm	19 - 22 mm	<b>1.19 mm</b>	<b>219.969.119</b>	<b>219.969.119S</b>

## Saw Blade fits to Aesculap® Acculan 3Ti

**90 MM LENGTH      12 - 19 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
90 mm	12 - 19 mm	<b>0.90 mm</b>	<b>219.999.090</b>	<b>219.999.090S</b>
90 mm	12 - 19 mm	<b>1.00 mm</b>	<b>219.999.100</b>	<b>219.999.100S</b>
90 mm	12 - 19 mm	<b>1.19 mm</b>	<b>219.999.119</b>	<b>219.999.119S</b>
90 mm	12 - 19 mm	<b>1.23 mm</b>	<b>219.999.123</b>	<b>219.999.123S</b>
90 mm	12 - 19 mm	<b>1.25 mm</b>	<b>219.999.125</b>	<b>219.999.125S</b>
90 mm	12 - 19 mm	<b>1.27 mm</b>	<b>219.999.127</b>	<b>219.999.127S</b>
90 mm	12 - 19 mm	<b>1.37 mm</b>	<b>219.999.137</b>	<b>219.999.137S</b>
90 mm	12 - 19 mm	<b>1.47 mm</b>	<b>219.999.147</b>	<b>219.999.147S</b>

**70 MM LENGTH      25 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
70 mm	25 mm	<b>0.90 mm</b>	<b>219.257.090</b>	<b>219.257.090S</b>
70 mm	25 mm	<b>1.00 mm</b>	<b>219.257.100</b>	<b>219.257.100S</b>
70 mm	25 mm	<b>1.19 mm</b>	<b>219.257.119</b>	<b>219.257.119S</b>
70 mm	25 mm	<b>1.23 mm</b>	<b>219.257.123</b>	<b>219.257.123S</b>
70 mm	25 mm	<b>1.25 mm</b>	<b>219.257.125</b>	<b>219.257.125S</b>
70 mm	25 mm	<b>1.27 mm</b>	<b>219.257.127</b>	<b>219.257.127S</b>
70 mm	25 mm	<b>1.37 mm</b>	<b>219.257.137</b>	<b>219.257.137S</b>
70 mm	25 mm	<b>1.47 mm</b>	<b>219.257.147</b>	<b>219.257.147S</b>

## Saw Blade fits to Aesculap® Acculan 3Ti

**70 MM LENGTH      19 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
70 mm	19 mm	<b>0.90 mm</b>	<b>219.197.090</b>	<b>219.197.090S</b>
70 mm	19 mm	<b>1.00 mm</b>	<b>219.197.100</b>	<b>219.197.100S</b>
70 mm	19 mm	<b>1.19 mm</b>	<b>219.197.119</b>	<b>219.197.119S</b>
70 mm	19 mm	<b>1.23 mm</b>	<b>219.197.123</b>	<b>219.197.123S</b>
70 mm	19 mm	<b>1.25 mm</b>	<b>219.197.125</b>	<b>219.197.125S</b>
70 mm	19 mm	<b>1.27 mm</b>	<b>219.197.127</b>	<b>219.197.127S</b>
70 mm	19 mm	<b>1.37 mm</b>	<b>219.197.137</b>	<b>219.197.137S</b>
70 mm	19 mm	<b>1.47 mm</b>	<b>219.197.147</b>	<b>219.197.147S</b>

**70 MM LENGTH      12 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
70 mm	12 mm	<b>1.27 mm</b>	<b>219.127.127</b>	<b>219.127.127S</b>

**50 MM LENGTH      12 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
50 mm	12 mm	<b>1.27 mm</b>	<b>219.125.127</b>	<b>219.125.127S</b>

## Saw Blade fits to MicroAire® (Series 7000)

**90 MM LENGTH      25 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
90 mm	25 mm	<b>0.90 mm</b>	<b>220.259.090</b>	<b>220.259.090S</b>
90 mm	25 mm	<b>1.00 mm</b>	<b>220.259.100</b>	<b>220.259.100S</b>
90 mm	25 mm	<b>1.19 mm</b>	<b>220.259.119</b>	<b>220.259.119S</b>
90 mm	25 mm	<b>1.23 mm</b>	<b>220.259.123</b>	<b>220.259.123S</b>
90 mm	25 mm	<b>1.25 mm</b>	<b>220.259.125</b>	<b>220.259.125S</b>
90 mm	25 mm	<b>1.27 mm</b>	<b>220.259.127</b>	<b>220.259.127S</b>
90 mm	25 mm	<b>1.37 mm</b>	<b>220.259.137</b>	<b>220.259.137S</b>
90 mm	25 mm	<b>1.47 mm</b>	<b>220.259.147</b>	<b>220.259.147S</b>

**90 MM LENGTH      19 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
90 mm	19 mm	<b>0.90 mm</b>	<b>220.199.090</b>	<b>220.199.090S</b>
90 mm	19 mm	<b>1.00 mm</b>	<b>220.199.100</b>	<b>220.199.100S</b>
90 mm	19 mm	<b>1.19 mm</b>	<b>220.199.119</b>	<b>220.199.119S</b>
90 mm	19 mm	<b>1.23 mm</b>	<b>220.199.123</b>	<b>220.199.123S</b>
90 mm	19 mm	<b>1.25 mm</b>	<b>220.199.125</b>	<b>220.199.125S</b>
90 mm	19 mm	<b>1.27 mm</b>	<b>220.199.127</b>	<b>220.199.127S</b>
90 mm	19 mm	<b>1.27 mm</b>	<b>220.199.127</b>	<b>220.199.127S</b>
90 mm	19 mm	<b>1.37 mm</b>	<b>220.199.137</b>	<b>220.199.137S</b>

## Saw Blade fits to MicroAire® (Series 7000)

90 MM LENGTH

12 MM WIDTH



Length	Width	Thickness	REF No non-sterile	REF No sterile
90 mm	12 mm	<b>0.90 mm</b>	<b>220.129.090</b>	<b>220.129.090S</b>
90 mm	12 mm	<b>1.00 mm</b>	<b>220.129.100</b>	<b>220.129.100S</b>
90 mm	12 mm	<b>1.19 mm</b>	<b>220.129.119</b>	<b>220.129.119S</b>
90 mm	12 mm	<b>1.23 mm</b>	<b>220.129.123</b>	<b>220.129.123S</b>
90 mm	12 mm	<b>1.25 mm</b>	<b>220.129.125</b>	<b>220.129.125S</b>
90 mm	12 mm	<b>1.27 mm</b>	<b>220.129.127</b>	<b>220.129.127S</b>
90 mm	12 mm	<b>1.37 mm</b>	<b>220.129.137</b>	<b>220.129.137S</b>
90 mm	12 mm	<b>1.47 mm</b>	<b>220.129.147</b>	<b>220.129.147S</b>

## Saw Blade fits to Arthrex Power System 600

**110 MM LENGTH    19 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
110 mm	19 mm	1.27 mm	221.191.127	221.191.127S

**90 MM LENGTH    19 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
90 mm	19 mm	1.27 mm	221.199.127	221.199.127S

**90 MM LENGTH    12 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
90 mm	12 mm	1.27 mm	221.129.127	221.129.127S

## Performance Line

Manufacturer: **GOMINA**

Cutting Thickness

Cutting length in mm




Swiss Quality

Reference-No


LOT-Nr.

Gomina Teeth-  
Design Performance Line


## HUB OVERVIEW




Aesculap®




Stryker®, System EHD,  
2000 / 4 / 5 / 6 / 7




Linvatec Hall®




Aesculap®  
Acculan 3Ti




AO (Synthes®,  
Sodem®)



MicroAire®  
(Series 7000)



deSoutter  
(Multidrive™)



Arthrex Power  
System 600

## Saw Blade fits to Aesculap®

**90 MM LENGTH      25 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
90 mm	25 mm	<b>0.90 mm</b>	<b>262.259.090</b>	<b>262.259.090S</b>
90 mm	25 mm	<b>1.00 mm</b>	<b>262.259.100</b>	<b>262.259.100S</b>
90 mm	25 mm	<b>1.19 mm</b>	<b>262.259.119</b>	<b>262.259.119S</b>
90 mm	25 mm	<b>1.27 mm</b>	<b>262.259.127</b>	<b>262.259.127S</b>
90 mm	25 mm	<b>1.37 mm</b>	<b>262.259.137</b>	<b>262.259.137S</b>
90 mm	25 mm	<b>1.47 mm</b>	<b>262.259.147</b>	<b>262.259.147S</b>

**90 MM LENGTH      19 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
90 mm	19 mm	<b>0.90 mm</b>	<b>262.199.090</b>	<b>262.199.090S</b>
90 mm	19 mm	<b>1.00 mm</b>	<b>262.199.100</b>	<b>262.199.100S</b>
90 mm	19 mm	<b>1.19 mm</b>	<b>262.199.119</b>	<b>262.199.119S</b>
90 mm	19 mm	<b>1.27 mm</b>	<b>262.199.127</b>	<b>262.199.127S</b>
90 mm	19 mm	<b>1.37 mm</b>	<b>262.199.137</b>	<b>262.199.137S</b>
90 mm	19 mm	<b>1.47 mm</b>	<b>262.199.147</b>	<b>262.199.147S</b>



## Saw Blade fits to Aesculap®

**90 MM LENGTH      12 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
90 mm	12 mm	<b>0.90 mm</b>	<b>262.129.090</b>	<b>262.129.090S</b>
90 mm	12 mm	<b>1.00 mm</b>	<b>262.129.100</b>	<b>262.129.100S</b>
90 mm	12 mm	<b>1.19 mm</b>	<b>262.129.119</b>	<b>262.129.119S</b>
90 mm	12 mm	<b>1.27 mm</b>	<b>262.129.127</b>	<b>262.129.127S</b>
90 mm	12 mm	<b>1.37 mm</b>	<b>262.129.137</b>	<b>262.129.137S</b>
90 mm	12 mm	<b>1.47 mm</b>	<b>262.129.147</b>	<b>262.129.147S</b>


**65 MM LENGTH      10 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
65 mm	10 mm	<b>1.27 mm</b>	<b>262.107.127</b>	<b>262.107.127S</b>

## Saw Blade fits to Linvatec Hall®

**110 MM LENGTH    19 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
110 mm	19 mm	<b>1.19 mm</b>	<b>263.191.119</b>	<b>263.191.119S</b>
110 mm	19 mm	<b>1.27 mm</b>	<b>263.191.127</b>	<b>263.191.127S</b>

**90 MM LENGTH    25 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
90 mm	25 mm	<b>0.90 mm</b>	<b>263.259.090</b>	<b>263.259.090S</b>
90 mm	25 mm	<b>1.00 mm</b>	<b>263.259.100</b>	<b>263.259.100S</b>
90 mm	25 mm	<b>1.19 mm</b>	<b>263.259.119</b>	<b>263.259.119S</b>
90 mm	25 mm	<b>1.27 mm</b>	<b>263.259.127</b>	<b>263.259.127S</b>
90 mm	25 mm	<b>1.37 mm</b>	<b>263.259.137</b>	<b>263.259.137S</b>
90 mm	25 mm	<b>1.47 mm</b>	<b>263.259.147</b>	<b>263.259.147S</b>

**90 MM LENGTH    19 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
90 mm	19 mm	<b>0.90 mm</b>	<b>263.199.090</b>	<b>263.199.090S</b>
90 mm	19 mm	<b>1.00 mm</b>	<b>263.199.100</b>	<b>263.199.100S</b>
90 mm	19 mm	<b>1.19 mm</b>	<b>263.199.119</b>	<b>263.199.119S</b>
90 mm	19 mm	<b>1.27 mm</b>	<b>263.199.127</b>	<b>263.199.127S</b>
90 mm	19 mm	<b>1.37 mm</b>	<b>263.199.137</b>	<b>263.199.137S</b>
90 mm	19 mm	<b>1.47 mm</b>	<b>263.199.147</b>	<b>263.199.147S</b>

## Saw Blade fits to Linvatec Hall®

**90 MM LENGTH      12 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
90 mm	12 mm	<b>0.90 mm</b>	<b>263.129.090</b>	<b>263.129.090S</b>
90 mm	12 mm	<b>1.00 mm</b>	<b>263.129.100</b>	<b>263.129.100S</b>
90 mm	12 mm	<b>1.19 mm</b>	<b>263.129.119</b>	<b>263.129.119S</b>
90 mm	12 mm	<b>1.27 mm</b>	<b>263.129.127</b>	<b>263.129.127S</b>
90 mm	12 mm	<b>1.37 mm</b>	<b>263.129.137</b>	<b>263.129.137S</b>
90 mm	12 mm	<b>1.47 mm</b>	<b>263.129.147</b>	<b>263.129.147S</b>

**70 MM LENGTH      25 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
70 mm	25 mm	<b>1.00 mm</b>	<b>263.257.100</b>	<b>263.257.100S</b>

**70 MM LENGTH      19 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
70 mm	19 mm	<b>0.90 mm</b>	<b>263.197.090</b>	<b>263.197.090S</b>

## Saw Blade fits to Linvatec Hall®

**70 MM LENGTH      12 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
70 mm	12 mm	<b>0.90 mm</b>	<b>263.127.090</b>	<b>263.127.090S</b>
70 mm	12 mm	<b>1.00 mm</b>	<b>263.127.100</b>	<b>263.127.100S</b>
70 mm	12 mm	<b>1.19 mm</b>	<b>263.127.119</b>	<b>263.127.119S</b>
70 mm	12 mm	<b>1.27 mm</b>	<b>263.127.127</b>	<b>263.127.127S</b>
70 mm	12 mm	<b>1.37 mm</b>	<b>263.127.137</b>	<b>263.127.137S</b>
70 mm	12 mm	<b>1.47 mm</b>	<b>263.127.147</b>	<b>263.127.147S</b>


**55 MM LENGTH      10 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
55 mm	10 mm	<b>1.27 mm</b>	<b>263.106.127</b>	<b>263.106.127S</b>


## Saw Blade fits to AO (Synthes®, Sodem®)

**110 MM LENGTH      19 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
110 mm	19 mm	<b>1.19 mm</b>	<b>265.191.119</b>	<b>265.191.119S</b>
110 mm	19 mm	<b>1.27 mm</b>	<b>265.191.127</b>	<b>265.191.127S</b>

**90 MM LENGTH      25 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
90 mm	25 mm	<b>0.90 mm</b>	<b>265.259.090</b>	<b>265.259.090S</b>
90 mm	25 mm	<b>1.00 mm</b>	<b>265.259.100</b>	<b>265.259.100S</b>
90 mm	25 mm	<b>1.19 mm</b>	<b>265.259.119</b>	<b>265.259.119S</b>
90 mm	25 mm	<b>1.27 mm</b>	<b>265.259.127</b>	<b>265.259.127S</b>
90 mm	25 mm	<b>1.37 mm</b>	<b>265.259.137</b>	<b>265.259.137S</b>
90 mm	25 mm	<b>1.47 mm</b>	<b>265.259.147</b>	<b>265.259.147S</b>

**90 MM LENGTH      19 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
90 mm	19 mm	<b>0.90 mm</b>	<b>265.199.090</b>	<b>265.199.090S</b>
90 mm	19 mm	<b>1.00 mm</b>	<b>265.199.100</b>	<b>265.199.100S</b>
90 mm	19 mm	<b>1.19 mm</b>	<b>265.199.119</b>	<b>265.199.119S</b>
90 mm	19 mm	<b>1.23 mm</b>	<b>265.199.123</b>	<b>265.199.123S</b>
90 mm	19 mm	<b>1.27 mm</b>	<b>265.199.127</b>	<b>265.199.127S</b>
90 mm	19 mm	<b>1.37 mm</b>	<b>265.199.137</b>	<b>265.199.137S</b>
90 mm	19 mm	<b>1.47 mm</b>	<b>265.199.147</b>	<b>265.199.147S</b>

## Saw Blade fits to AO (Synthes®, Sodem®)

**90 MM LENGTH      12 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
90 mm	12 mm	<b>0.90 mm</b>	<b>265.129.090</b>	<b>265.129.090S</b>
90 mm	12 mm	<b>1.00 mm</b>	<b>265.129.100</b>	<b>265.129.100S</b>
90 mm	12 mm	<b>1.19 mm</b>	<b>265.129.119</b>	<b>265.129.119S</b>
90 mm	12 mm	<b>1.23 mm</b>	<b>265.129.123</b>	<b>265.129.123S</b>
90 mm	12 mm	<b>1.27 mm</b>	<b>265.129.127</b>	<b>265.129.127S</b>
90 mm	12 mm	<b>1.37 mm</b>	<b>265.129.137</b>	<b>265.129.137S</b>
90 mm	12 mm	<b>1.47 mm</b>	<b>265.129.147</b>	<b>265.129.147S</b>

**80 MM LENGTH      12 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
80 mm	12 mm	<b>1.27 mm</b>	<b>265.128.127</b>	<b>265.128.127S</b>

**80 MM LENGTH      8 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
80 mm	8 mm	<b>1.27 mm</b>	<b>265.088.127</b>	<b>265.088.127S</b>

## Saw Blade fits to AO (Synthes®, Sodem®)

**70 MM LENGTH      19 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
70 mm	19 mm	<b>0.90 mm</b>	<b>265.197.090</b>	<b>265.197.090S</b>
70 mm	19 mm	<b>1.00 mm</b>	<b>265.197.100</b>	<b>265.197.100S</b>
70 mm	19 mm	<b>1.19 mm</b>	<b>265.197.119</b>	<b>265.197.119S</b>
70 mm	19 mm	<b>1.27 mm</b>	<b>265.197.127</b>	<b>265.197.127S</b>
70 mm	19 mm	<b>1.37 mm</b>	<b>265.197.137</b>	<b>265.197.137S</b>
70 mm	19 mm	<b>1.47 mm</b>	<b>265.197.147</b>	<b>265.197.147S</b>

**70 MM LENGTH      12 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
70 mm	12 mm	<b>0.90 mm</b>	<b>265.127.090</b>	<b>265.127.090S</b>
70 mm	12 mm	<b>1.00 mm</b>	<b>265.127.100</b>	<b>265.127.100S</b>
70 mm	12 mm	<b>1.19 mm</b>	<b>265.127.119</b>	<b>265.127.119S</b>
70 mm	12 mm	<b>1.27 mm</b>	<b>265.127.127</b>	<b>265.127.127S</b>
70 mm	12 mm	<b>1.37 mm</b>	<b>265.127.137</b>	<b>265.127.137S</b>
70 mm	12 mm	<b>1.47 mm</b>	<b>265.127.147</b>	<b>265.127.147S</b>

**70 MM LENGTH      8 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
70 mm	8 mm	<b>1.27 mm</b>	<b>265.087.127</b>	<b>265.087.127S</b>

## Saw Blade fits to AO (Synthes®, Sodem®)

**65 MM LENGTH      10 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
65 mm	10 mm	1.27 mm	265.107.127	265.107.127S



## Saw Blade fits to deSoutter (Multidrive™)

**100 MM LENGTH    25 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
100 mm	25 mm	1.27 mm	267.250.127	267.250.127S

**90 MM LENGTH    25 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
90 mm	25 mm	0.90 mm	267.259.090	267.259.090S
90 mm	25 mm	1.00 mm	267.259.100	267.259.100S
90 mm	25 mm	1.19 mm	267.259.119	267.259.119S
90 mm	25 mm	1.27 mm	267.259.127	267.259.127S
90 mm	25 mm	1.37 mm	267.259.137	267.259.137S
90 mm	25 mm	1.47 mm	267.259.147	267.259.147S

**90 MM LENGTH    19 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
90 mm	19 mm	0.90 mm	267.199.090	267.199.090S
90 mm	19 mm	1.00 mm	267.199.100	267.199.100S
90 mm	19 mm	1.19 mm	267.199.119	267.199.119S
90 mm	19 mm	1.27 mm	267.199.127	267.199.127S
90 mm	19 mm	1.37 mm	267.199.137	267.199.137S
90 mm	19 mm	1.47 mm	267.199.147	267.199.147S

## Saw Blade fits to deSoutter (Multidrive™)

**90 MM LENGTH      12 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
90 mm	12 mm	<b>0.90 mm</b>	<b>267.129.090</b>	<b>267.129.090S</b>
90 mm	12 mm	<b>1.00 mm</b>	<b>267.129.100</b>	<b>267.129.100S</b>
90 mm	12 mm	<b>1.19 mm</b>	<b>267.129.119</b>	<b>267.129.119S</b>
90 mm	12 mm	<b>1.27 mm</b>	<b>267.129.127</b>	<b>267.129.127S</b>
90 mm	12 mm	<b>1.37 mm</b>	<b>267.129.137</b>	<b>267.129.137S</b>
90 mm	12 mm	<b>1.47 mm</b>	<b>267.129.147</b>	<b>267.129.147S</b>

**80 MM LENGTH      12 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
80 mm	12 mm	<b>1.27 mm</b>	<b>267.128.127</b>	<b>267.128.127S</b>

**80 MM LENGTH      8 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
80 mm	8 mm	<b>1.27 mm</b>	<b>267.088.127</b>	<b>267.088.127S</b>

## Saw Blade fits to deSoutter (Multidrive™)

**70 MM LENGTH      25 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
70 mm	25 mm	<b>0.90 mm</b>	<b>267.257.090</b>	<b>267.257.090S</b>

**50 MM LENGTH      12 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
50 mm	12 mm	<b>0.80 mm</b>	<b>267.125.080</b>	<b>267.125.080S</b>

## Saw Blade fits to Stryker®, System EHD, 2000 / 4 / 5 / 6

**110 MM LENGTH    25 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
110 mm	25 mm	1.27 mm	268.251.127	268.251.127S

**90 MM LENGTH    25 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
90 mm	25 mm	0.90 mm	268.259.090	268.259.090S
90 mm	25 mm	1.00 mm	268.259.100	268.259.100S
90 mm	25 mm	1.19 mm	268.259.119	268.259.119S
90 mm	25 mm	1.27 mm	268.259.127	268.259.127S
90 mm	25 mm	1.37 mm	268.259.137	268.259.137S
90 mm	25 mm	1.47 mm	268.259.147	268.259.147S


**90 MM LENGTH    19 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
90 mm	19 mm	0.90 mm	268.199.090	268.199.090S
90 mm	19 mm	1.00 mm	268.199.100	268.199.100S
90 mm	19 mm	1.19 mm	268.199.119	268.199.119S
90 mm	19 mm	1.23 mm	268.199.123	268.199.123S
90 mm	19 mm	1.25 mm	268.199.125	268.199.125S
90 mm	19 mm	1.27 mm	268.199.127	268.199.127S
90 mm	19 mm	1.37 mm	268.199.137	268.199.137S
90 mm	19 mm	1.47 mm	268.199.147	268.199.147S

## Saw Blade fits to Stryker®, System EHD, 2000 / 4 / 5 / 6

**90 MM LENGTH      12 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
90 mm	12 mm	<b>0.90 mm</b>	<b>268.129.090</b>	<b>268.129.090S</b>
90 mm	12 mm	<b>1.00 mm</b>	<b>268.129.100</b>	<b>268.129.100S</b>
90 mm	12 mm	<b>1.19 mm</b>	<b>268.129.119</b>	<b>268.129.119S</b>
90 mm	12 mm	<b>1.27 mm</b>	<b>268.129.127</b>	<b>268.129.127S</b>
90 mm	12 mm	<b>1.37 mm</b>	<b>268.129.137</b>	<b>268.129.137S</b>
90 mm	12 mm	<b>1.47 mm</b>	<b>268.129.147</b>	<b>268.129.147S</b>

**80 MM LENGTH      12 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
80 mm	12 mm	<b>1.27 mm</b>	<b>268.128.127</b>	<b>268.128.127S</b>

**80 MM LENGTH      8 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
80 mm	8 mm	<b>1.27 mm</b>	<b>268.088.127</b>	<b>268.088.127S</b>

## Saw Blade fits to Stryker®, System EHD, 2000 / 4 / 5 / 6

**70 MM LENGTH      25 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
70 mm	25 mm	<b>0.90 mm</b>	<b>268.257.090</b>	<b>268.257.090S</b>
70 mm	25 mm	<b>1.00 mm</b>	<b>268.257.100</b>	<b>268.257.100S</b>
70 mm	25 mm	<b>1.19 mm</b>	<b>268.257.119</b>	<b>268.257.119S</b>
70 mm	25 mm	<b>1.27 mm</b>	<b>268.257.127</b>	<b>268.257.127S</b>
70 mm	25 mm	<b>1.37 mm</b>	<b>268.257.137</b>	<b>268.257.137S</b>
70 mm	25 mm	<b>1.47 mm</b>	<b>268.257.147</b>	<b>268.257.147S</b>

**70 MM LENGTH      19 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
70 mm	19 mm	<b>0.90 mm</b>	<b>268.197.090</b>	<b>268.197.090S</b>

**70 MM LENGTH      12 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
70 mm	12 mm	<b>0.90 mm</b>	<b>268.127.090</b>	<b>268.127.090S</b>
70 mm	12 mm	<b>1.00 mm</b>	<b>268.127.100</b>	<b>268.127.100S</b>
70 mm	12 mm	<b>1.19 mm</b>	<b>268.127.119</b>	<b>268.127.119S</b>
70 mm	12 mm	<b>1.27 mm</b>	<b>268.127.127</b>	<b>268.127.127S</b>
70 mm	12 mm	<b>1.37 mm</b>	<b>268.127.137</b>	<b>268.127.137S</b>
70 mm	12 mm	<b>1.47 mm</b>	<b>268.127.147</b>	<b>268.127.147S</b>

**Saw Blade fits to Stryker®, System EHD, 2000 / 4 / 5 / 6**

**55 MM LENGTH      10 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
55 mm	10 mm	1.27 mm	268.106.127	268.106.127S

## Saw Blade fits to Aesculap® Acculan 3Ti

**90 MM LENGTH      25 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
90 mm	25 mm	<b>0.90 mm</b>	<b>269.259.090</b>	<b>269.259.090S</b>
90 mm	25 mm	<b>1.00 mm</b>	<b>269.259.100</b>	<b>269.259.100S</b>
90 mm	25 mm	<b>1.19 mm</b>	<b>269.259.119</b>	<b>269.259.119S</b>
90 mm	25 mm	<b>1.27 mm</b>	<b>269.259.127</b>	<b>269.259.127S</b>
90 mm	25 mm	<b>1.37 mm</b>	<b>269.259.137</b>	<b>269.259.137S</b>
90 mm	25 mm	<b>1.47 mm</b>	<b>269.259.147</b>	<b>269.259.147S</b>

**90 MM LENGTH      19 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
90 mm	19 mm	<b>0.90 mm</b>	<b>269.199.090</b>	<b>269.199.090S</b>
90 mm	19 mm	<b>1.00 mm</b>	<b>269.199.100</b>	<b>269.199.100S</b>
90 mm	19 mm	<b>1.19 mm</b>	<b>269.199.119</b>	<b>269.199.119S</b>
90 mm	19 mm	<b>1.27 mm</b>	<b>269.199.127</b>	<b>269.199.127S</b>
90 mm	19 mm	<b>1.37 mm</b>	<b>269.199.137</b>	<b>269.199.137S</b>
90 mm	19 mm	<b>1.47 mm</b>	<b>269.199.147</b>	<b>269.199.147S</b>



## Saw Blade fits to Aesculap® Acculan 3Ti

**90 MM LENGTH      12 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
90 mm	12 mm	<b>0.90 mm</b>	<b>269.129.090</b>	<b>269.129.090S</b>
90 mm	12 mm	<b>1.00 mm</b>	<b>269.129.100</b>	<b>269.129.100S</b>
90 mm	12 mm	<b>1.19 mm</b>	<b>269.129.119</b>	<b>269.129.119S</b>
90 mm	12 mm	<b>1.27 mm</b>	<b>269.129.127</b>	<b>269.129.127S</b>
90 mm	12 mm	<b>1.37 mm</b>	<b>269.129.137</b>	<b>269.129.137S</b>
90 mm	12 mm	<b>1.47 mm</b>	<b>269.129.147</b>	<b>269.129.147S</b>

**80 MM LENGTH      12 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
80 mm	12 mm	<b>1.27 mm</b>	<b>269.128.127</b>	<b>269.128.127S</b>

**80 MM LENGTH      8 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
80 mm	8 mm	<b>1.27 mm</b>	<b>269.088.127</b>	<b>269.088.127S</b>

## Saw Blade fits to Aesculap® Acculan 3Ti

**65 MM LENGTH      10 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
65 mm	10 mm	1.27 mm	269.107.127	269.107.127S

## Saw Blade fits to MicroAire® (Series 7000)

**110 MM LENGTH    19 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
110 mm	19 mm	1.19 mm	270.191.119	270.191.119S

**90 MM LENGTH    25 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
90 mm	25 mm	1.27 mm	270.259.127	270.259.127S

**90 MM LENGTH    19 MM WIDTH**



Length	Width	Thickness	REF No non-sterile	REF No sterile
90 mm	19 mm	0.90 mm	270.199.090	270.199.090S
90 mm	19 mm	1.00 mm	270.199.100	270.199.100S
90 mm	19 mm	1.19 mm	270.199.119	270.199.119S
90 mm	19 mm	1.27 mm	270.199.127	270.199.127S

## Saw Blade fits to MicroAire® (Series 7000)

90 MM LENGTH


12 MM WIDTH




Length	Width	Thickness	REF No non-sterile	REF No sterile
90 mm	12 mm	1.00 mm	270.129.100	270.129.100S
90 mm	12 mm	1.19 mm	270.129.119	270.129.119S
90 mm	12 mm	1.27 mm	270.129.127	270.129.127S

## Reciprocating Saw Blades


### HUB OVERVIEW



Linvatec Hall®,  
Power Pro,  
Versipower




Synthes®




Stryker®,  
System EHD,  
2000 / 4 / 5 / 6 / 7



Zimmer®



Tuke®



Aesculap®  
Acculan 3Ti

## Reciprocating Saw Blade fits to Linvatec Hall®, Power Pro, Versipower



Length	Width	Thickness	REF No non-sterile	REF No sterile
70 mm	12 mm	1.10 mm	231.270.110	231.270.110S



Length	Width	Thickness	REF No non-sterile	REF No sterile
60 mm	9 mm	1.00 mm	231.160.100	231.160.100S



Length	Width	Thickness	REF No non-sterile	REF No sterile
77 mm	12 mm	1.10 mm	231.177.110	231.177.110S



Length	Width	Thickness	REF No non-sterile	REF No sterile
85 mm	12 mm	1.20 mm	231.785.120	231.785.120S

## STERNUM



Length	Width	Thickness	REF No non-sterile	REF No sterile
34 mm	9.50 mm	1.00 mm	231.334.100	231.334.100S

## Reciprocating Saw Blade fits to Stryker®, System EHD, 2000 / 4 / 5 / 6 / 7



Length	Width	Thickness	REF No non-sterile	REF No sterile
70 mm	12 mm	1.10 mm	232.270.110	232.270.110S



Length	Width	Thickness	REF No non-sterile	REF No sterile
60 mm	9 mm	1.00 mm	232.160.100	232.160.100S



Length	Width	Thickness	REF No non-sterile	REF No sterile
77 mm	12 mm	1.10 mm	232.177.110	232.177.110S
77 mm	12 mm	1.23 mm	232.177.123	232.177.123S



Length	Width	Thickness	REF No non-sterile	REF No sterile
55 mm	10 mm	0.85 mm	232.755.085	232.755.085S



Length	Width	Thickness	REF No non-sterile	REF No sterile
85 mm	12 mm	1.20 mm	232.785.120	232.785.120S

## Reciprocating Saw Blade fits to Tuke®




Length	Width	Thickness	REF No non-sterile	REF No sterile
65 mm	12.35 mm	0.80 mm	236.265.080	236.265.080S




Length	Width	Thickness	REF No non-sterile	REF No sterile
65 mm	12.35 mm	1.00 mm	236.165.100	236.165.100S



Length	Width	Thickness	REF No non-sterile	REF No sterile
66 mm	16.6 mm	0.80 mm	236.266.080	236.266.080S



Length	Width	Thickness	REF No non-sterile	REF No sterile
80 mm	12.35 mm	1.00 mm	236.180.100	236.180.100S



Length	Width	Thickness	REF No non-sterile	REF No sterile
87 mm	12.35 mm	0.80 mm	236.287.080	236.287.080S



## Reciprocating Saw Blade fits to Synthes®




Length	Width	Thickness	REF No non-sterile	REF No sterile
70 mm	12 mm	1.10 mm	235.270.110	235.270.110S



Length	Width	Thickness	REF No non-sterile	REF No sterile
60 mm	9 mm	1.05 mm	235.160.105	235.160.105S
60 mm	9 mm	1.10 mm	235.160.110	235.160.110S



Length	Width	Thickness	REF No non-sterile	REF No sterile
70 mm	10 mm	1.30 mm	235.170.130	235.170.130S



Length	Width	Thickness	REF No non-sterile	REF No sterile
77 mm	12 mm	1.10 mm	235.177.110	235.177.110S




Length	Width	Thickness	REF No non-sterile	REF No sterile
55 mm	10 mm	0.85 mm	235.755.085	235.755.085S



Length	Width	Thickness	REF No non-sterile	REF No sterile
85 mm	12 mm	1.20 mm	235.785.120	235.785.120S


11/2017

## Reciprocating Saw Blade fits to Synthes®



Length	Width	Thickness	REF No non-sterile	REF No sterile
80 mm	10 mm	1.40 mm	235.980.990	235.980.990S

## Reciprocating Saw Blade fits to Synthes®




Length	Width	Thickness	REF No non-sterile	REF No sterile
82 mm	10 mm	1.00 mm	245.182.100	245.182.100S



Length	Width	Thickness	REF No non-sterile	REF No sterile
86 mm	12 mm	1.20 mm	245.186.120	245.186.120S



Length	Width	Thickness	REF No non-sterile	REF No sterile
65 mm	12 mm	1.20 mm	245.165.120	245.165.120S



Length	Width	Thickness	REF No non-sterile	REF No sterile
88 mm	14 mm	1.20 mm	245.188.120	245.188.120S

## Reciprocating Saw Blade fits to Zimmer®



Length	Width	Thickness	REF No non-sterile	REF No sterile
73 mm	13 mm	0.90 mm	234.873.090	234.873.090S



Length	Width	Thickness	REF No non-sterile	REF No sterile
85 mm	12 mm	1.20 mm	234.785.120	234.785.120S


## Reciprocating Saw Blade fits to Aesculap® Acculan 3Ti



Length	Width	Thickness	REF No non-sterile	REF No sterile
70 mm	12 mm	1.10 mm	233.270.110	233.270.110S



Length	Width	Thickness	REF No non-sterile	REF No sterile
60 mm	9 mm	1.00 mm	233.160.100	233.160.100S



Length	Width	Thickness	REF No non-sterile	REF No sterile
77 mm	12 mm	1.10 mm	233.177.110	233.177.110S





Length	Width	Thickness	REF No non-sterile	REF No sterile
82 mm	10 mm	1.00 mm	233.782.100	233.782.100S




Length	Width	Thickness	REF No non-sterile	REF No sterile
85 mm	12 mm	1.20 mm	233.785.120	233.785.120S

## Reciprocating Keel Blade

### HUB OVERVIEW



## Reciprocating Keel Blade fits to Linvatec Hall®, Power Pro, Versipower



Reciprocating Keel Blade	REF No non-sterile	REF No sterile
cementless	not available	231.018.275S
cemented	not available	231.018.290S

## Keel Blade Set fits to Linvatec Hall®, Power Pro, Versipower


Keel Blade
cementless
cemented

Length	Width	Thickness
90 mm	12 mm	<b>0.90 mm</b>

Length	Width	Thickness
85 mm	12 mm	<b>1.20 mm</b>


Keel Blade Set	REF No non-sterile	REF No sterile
Keel Blade Set cementless	not available	231.000.001S
Keel Blade Set cemented	not available	231.000.002S

## Reciprocating Keel Blade fits to Stryker®, System EHD, 2000 / 4 / 5 / 6 / 7



Reciprocating Keel Blade	REF No non-sterile	REF No sterile
cementless	not available	232.018.275S
cemented	not available	232.018.290S

## Keel Blade Set fits to Stryker®, System EHD, 2000 / 4 / 5 / 6 / 7



Keel Blade
cementless
cemented


Length	Width	Thickness
90 mm	12 mm	0.90 mm

Length	Width	Thickness
85 mm	12 mm	1.20 mm

Keel Blade Set	REF No non-sterile	REF No sterile
Keel Blade Set cementless	not available	232.000.001S
Keel Blade Set cemented	not available	232.000.002S



## Reciprocating Keel Blade fits to Synthes®



Reciprocating Keel Blade	REF No non-sterile	REF No sterile
cementless	not available	235.018.275S
cemented	not available	235.018.290S

## Keel Blade Set fits to AO (Synthes®)


Keel Blade		
cementless		
cemented		

Length	Width	Thickness
90 mm	12 mm	<b>0.90 mm</b>

Length	Width	Thickness
85 mm	12 mm	<b>1.20 mm</b>

Keel Blade Set	REF No non-sterile	REF No sterile
Keel Blade Set cementless	not available	235.000.001S
Keel Blade Set cemented	not available	235.000.002S

## Reciprocating Keel Blade fits to Aesculap® Acculan 3Ti



Reciprocating Keel Blade	REF No non-sterile	REF No sterile
cementless	not available	233.018.275S
cemented	not available	233.018.290S

## Keel Blade Set fits to Aesculap® Acculan 3Ti


Keel Blade
cementless
cemented

Length	Width	Thickness
90 mm	12 mm	0.90 mm

Length	Width	Thickness
85 mm	12 mm	1.20 mm

Keel Blade Set	REF No non-sterile	REF No sterile
Keel Blade Set cementless	not available	233.000.001S
Keel Blade Set cemented	not available	233.000.002S

## Reciprocating Keel Blade fits to Zimmer®



Reciprocating Keel Blade	REF No non-sterile	REF No sterile
cementless	not available	234.018.275S
cemented	not available	234.018.290S

## Keel Blade Set fits to Zimmer®

Keel Blade
cementless
cemented

Length	Width	Thickness
90 mm	12 mm	0.90 mm


  

Length	Width	Thickness
85 mm	12 mm	1.20 mm


Keel Blade Set	REF No non-sterile	REF No sterile
Keel Blade Set cementless	not available	234.000.001S
Keel Blade Set cemented	not available	234.000.002S

## Saw Blades for Trauma


### HUB OVERVIEW




Aesculap®




deSoutter  
(Multidrive™)




Linvatec Hall®



Stryker®, System EHD,  
2000 / 4 / 5 / 6 / 7



AO (Synthes®,  
Sodem®)



Aesculap®  
Acculan 3Ti

## Saw Blade fits to Aesculap®



Length	Width	Thickness	REF No non-sterile	REF No sterile
35 mm	10 mm	<b>0.60 mm</b>	<b>202.1035.060</b>	<b>202.1035.060S</b>



Length	Width	Thickness	REF No non-sterile	REF No sterile
35 mm	20 mm	<b>0.60 mm</b>	<b>202.2035.060</b>	<b>202.2035.060S</b>



Length	Width	Thickness	REF No non-sterile	REF No sterile
48 mm	14 mm	<b>0.60 mm</b>	<b>202.1448.060</b>	<b>202.1448.060S</b>
48 mm	14 mm	<b>0.80 mm</b>	<b>202.1448.080</b>	<b>202.1448.080S</b>

## Saw Blade fits to Aesculap®



Length	Width	Thickness	REF No non-sterile	REF No sterile
48 mm	20 mm	<b>0.60 mm</b>	<b>202.2048.060</b>	<b>202.2048.060S</b>
48 mm	20 mm	<b>0.80 mm</b>	<b>202.2048.080</b>	<b>202.2048.080S</b>
48 mm	20 mm	<b>1.00 mm</b>	<b>202.2048.100</b>	<b>202.2048.100S</b>



Length	Width	Thickness	REF No non-sterile	REF No sterile
48 mm	27 mm	<b>0.60 mm</b>	<b>202.2748.060</b>	<b>202.2748.060S</b>
48 mm	27 mm	<b>0.80 mm</b>	<b>202.2748.080</b>	<b>202.2748.080S</b>



Length	Width	Thickness	REF No non-sterile	REF No sterile
70 mm	20 mm	<b>0.80 mm</b>	<b>202.2070.080</b>	<b>202.2070.080S</b>



Length	Width	Thickness	REF No non-sterile	REF No sterile
70 mm	27 mm	<b>0.80 mm</b>	<b>202.2770.080</b>	<b>202.2770.080S</b>

## Saw Blade fits to Linvatec Hall®



Length	Width	Thickness	REF No non-sterile	REF No sterile
63 mm	49 mm	0.90 mm	203.4963.090	203.4963.090S



Length	Width	Thickness	REF No non-sterile	REF No sterile
70 mm	20 mm	0.80 mm	203.2070.080	203.2070.080S




Length	Width	Thickness	REF No non-sterile	REF No sterile
70 mm	27 mm	0.80 mm	203.2770.080	203.2770.080S

## Saw Blade fits to AO (Synthes®, Sodem®)



Length	Width	Thickness	REF No non-sterile	REF No sterile
25 mm	6 mm	0.60 mm	206.0625.060	206.0625.060S



Length	Width	Thickness	REF No non-sterile	REF No sterile
25 mm	10 mm	0.60 mm	206.1025.060	206.1025.060S



Length	Width	Thickness	REF No non-sterile	REF No sterile
25 mm	14 mm	0.60 mm	206.1425.060	206.1425.060S



## Saw Blade fits to AO (Synthes®, Sodem®)



Length	Width	Thickness	REF No non-sterile	REF No sterile
48 mm	10 mm	0.60 mm	206.1048.060	206.1048.060S



Length	Width	Thickness	REF No non-sterile	REF No sterile
48 mm	14 mm	0.60 mm	206.1448.060	206.1448.060S



Length	Width	Thickness	REF No non-sterile	REF No sterile
48 mm	27 mm	0.60 mm	206.2748.060	206.2748.060S



Length	Width	Thickness	REF No non-sterile	REF No sterile
50 mm	20 mm	0.80 mm	206.2050.080	206.2050.080S

## Saw Blade fits to AO (Synthes®, Sodem®)



Length	Width	Thickness	REF No non-sterile	REF No sterile
70 mm	18 mm	<b>0.80 mm</b>	<b>206.1870.080</b>	<b>206.1870.080S</b>
70 mm	18 mm	<b>1.00 mm</b>	<b>206.1870.100</b>	<b>206.1870.100S</b>
70 mm	18 mm	<b>1.20 mm</b>	<b>206.1870.120</b>	<b>206.1870.120S</b>



Length	Width	Thickness	REF No non-sterile	REF No sterile
70 mm	18 mm	<b>1.20 mm</b>	<b>206.1870.120G</b>	<b>206.1870.120GS</b>



Length	Width	Thickness	REF No non-sterile	REF No sterile
70 mm	27 mm	<b>0.80 mm</b>	<b>206.2770.080</b>	<b>206.2770.080S</b>
70 mm	27 mm	<b>1.00 mm</b>	<b>206.2770.100</b>	<b>206.2770.100S</b>
70 mm	27 mm	<b>1.20 mm</b>	<b>206.2770.120</b>	<b>206.2770.120S</b>



Length	Width	Thickness	REF No non-sterile	REF No sterile
70 mm	27 mm	<b>1.20 mm</b>	<b>206.2770.120G</b>	<b>206.2770.120GS</b>

## Saw Blade fits to AO (Synthes®, Sodem®)



Length	Width	Thickness	REF No non-sterile	REF No sterile
90 mm	18 mm	<b>0.80 mm</b>	<b>206.1890.080</b>	<b>206.1890.080S</b>

## Saw Blade fits to deSoutter (Multidrive™)



Length	Width	Thickness	REF No non-sterile	REF No sterile
64 mm	34 mm	1.00 mm	207.3464.100	207.3464.100S

## Saw Blade fits to Stryker®, System EHD, 2000 / 4 / 5 / 6 / 7



Length	Width	Thickness	REF No non-sterile	REF No sterile
24 mm	9 mm	0.80 mm	208.0924.080	208.0924.080S



Length	Width	Thickness	REF No non-sterile	REF No sterile
25 mm	6 mm	0.60 mm	208.0625.060	208.0625.060S



Length	Width	Thickness	REF No non-sterile	REF No sterile
25 mm	14 mm	0.60 mm	208.1425.060	208.1425.060S



Length	Width	Thickness	REF No non-sterile	REF No sterile
35 mm	9 mm	0.80 mm	208.0935.080	208.0935.080S



Length	Width	Thickness	REF No non-sterile	REF No sterile
35 mm	20 mm	1.10 mm	208.2035.110	208.2035.110S

## Saw Blade fits to Stryker®, System EHD, 2000 / 4 / 5 / 6 / 7



Length	Width	Thickness	REF No non-sterile	REF No sterile
48 mm	14 mm	<b>0.60 mm</b>	<b>208.1448.060</b>	<b>208.1448.060S</b>



Length	Width	Thickness	REF No non-sterile	REF No sterile
48 mm	27 mm	<b>0.60 mm</b>	<b>208.2748.060</b>	<b>208.2748.060S</b>
48 mm	27 mm	<b>0.80 mm</b>	<b>208.2748.080</b>	<b>208.2748.080S</b>
48 mm	27 mm	<b>1.10 mm</b>	<b>208.2748.110</b>	<b>208.2748.110S</b>



Length	Width	Thickness	REF No non-sterile	REF No sterile
60 mm	30 mm	<b>0.80 mm</b>	<b>208.3060.080</b>	<b>208.3060.080S</b>

## Saw Blade fits to Stryker®, System EHD, 2000 / 4 / 5 / 6 / 7



Length	Width	Thickness	REF No non-sterile	REF No sterile
64 mm	34 mm	0.80 mm	208.3464.080	208.3464.080S




Length	Width	Thickness	REF No non-sterile	REF No sterile
70 mm	27 mm	0.80 mm	208.2770.080	208.2770.080S
70 mm	27 mm	1.20 mm	208.2770.120	208.2770.120S



Length	Width	Thickness	REF No non-sterile	REF No sterile
70 mm	18 mm	0.80 mm	208.1870.080	208.1870.080S


## Saw Blade fits to Aesculap® Acculan 3Ti




Length	Width	Thickness	REF No non-sterile	REF No sterile
25 mm	6 mm	0.60 mm	209.0625.060	209.0625.060S



Length	Width	Thickness	REF No non-sterile	REF No sterile
25 mm	14 mm	0.60 mm	209.1425.060	209.1425.060S



Length	Width	Thickness	REF No non-sterile	REF No sterile
35 mm	10 mm	0.60 mm	209.1035.060	209.1035.060S
35 mm	10 mm	0.80 mm	209.1035.080	209.1035.080S



Length	Width	Thickness	REF No non-sterile	REF No sterile
35 mm	14 mm	0.60 mm	209.1435.060	209.1435.060S



## Saw Blade fits to Aesculap® Acculan 3Ti



Length	Width	Thickness	REF No non-sterile	REF No sterile
35 mm	20 mm	<b>0.60 mm</b>	<b>209.2035.060</b>	<b>209.2035.060S</b>
35 mm	20 mm	<b>1.10 mm</b>	<b>209.2035.110</b>	<b>209.2035.110S</b>



Length	Width	Thickness	REF No non-sterile	REF No sterile
48 mm	14 mm	<b>0.60 mm</b>	<b>209.1448.060</b>	<b>209.1448.060S</b>
48 mm	14 mm	<b>1.00 mm</b>	<b>209.1448.100</b>	<b>209.1448.100S</b>



Length	Width	Thickness	REF No non-sterile	REF No sterile
48 mm	20 mm	<b>0.60 mm</b>	<b>209.2048.060</b>	<b>209.2048.060S</b>
48 mm	20 mm	<b>0.80 mm</b>	<b>209.2048.080</b>	<b>209.2048.080S</b>



Length	Width	Thickness	REF No non-sterile	REF No sterile
48 mm	27 mm	<b>0.60 mm</b>	<b>209.2748.060</b>	<b>209.2748.060S</b>
48 mm	27 mm	<b>0.80 mm</b>	<b>209.2748.080</b>	<b>209.2748.080S</b>
48 mm	27 mm	<b>1.10 mm</b>	<b>209.2748.110</b>	<b>209.2748.110S</b>

## Saw Blade fits to Aesculap® Acculan 3Ti



Length	Width	Thickness	REF No non-sterile	REF No sterile
70 mm	18 mm	<b>0.80 mm</b>	<b>209.1870.080</b>	<b>209.1870.080S</b>




Length	Width	Thickness	REF No non-sterile	REF No sterile
70 mm	20 mm	<b>0.80 mm</b>	<b>209.2070.080</b>	<b>209.2070.080S</b>
70 mm	20 mm	<b>1.10 mm</b>	<b>209.2070.110</b>	<b>209.2070.110S</b>




Length	Width	Thickness	REF No non-sterile	REF No sterile
70 mm	27 mm	<b>0.80 mm</b>	<b>209.2770.080</b>	<b>209.2770.080S</b>
70 mm	27 mm	<b>1.20 mm</b>	<b>209.2770.120</b>	<b>209.2770.120S</b>

## Saw Blades for small Bones


### HUB OVERVIEW




Stryker®, TPS




Colibri (Synthes®)




Smith & Nephew  
DYONICS™ Saw REF  
7205786 and 7205791



EPEN (Synthes®)




Linvatec Hall®,  
MicroChoice®,  
Micro100™, Micro E®



Arthrex



Linvatec Hall®,  
MicroChoice®,  
Micro100™, Micro E®



Aesculap® Acculan 3Ti

## Saw Blade fits to Stryker®, TPS



Length	Width	Thickness	REF No non-sterile	REF No sterile
18 mm	5 mm	0.40 mm	281.0518.040	281.0518.040S



Length	Width	Thickness	REF No non-sterile	REF No sterile
31 mm	5 mm	0.40 mm	281.0531.040	281.0531.040S



Length	Width	Thickness	REF No non-sterile	REF No sterile
15 mm	7 mm	0.40 mm	281.0715.040	281.0715.040S



Length	Width	Thickness	REF No non-sterile	REF No sterile
21 mm	7 mm	0.40 mm	281.0721.040	281.0721.040S



Length	Width	Thickness	REF No non-sterile	REF No sterile
31 mm	7 mm	0.40 mm	281.0731.040	281.0731.040S
31 mm	7 mm	0.60 mm	281.0731.060	281.0731.060S



Length	Width	Thickness	REF No non-sterile	REF No sterile
18 mm	9 mm	0.40 mm	281.0918.040	281.0918.040S



Length	Width	Thickness	REF No non-sterile	REF No sterile
25 mm	9 mm	0.40 mm	281.0925.040	281.0925.040S



Length	Width	Thickness	REF No non-sterile	REF No sterile
31 mm	9 mm	0.40 mm	281.0931.040	281.0931.040S
31 mm	9 mm	0.60 mm	281.0931.060	281.0931.060S



Length	Width	Thickness	REF No non-sterile	REF No sterile
35 mm	9 mm	0.50 mm	281.0935.050	281.0935.050S
35 mm	9 mm	0.60 mm	281.0935.060	281.0935.060S



Length	Width	Thickness	REF No non-sterile	REF No sterile
40 mm	11 mm	0.50 mm	281.1140.050	281.1140.050S

11/2017

## Saw Blade fits to Stryker®, TPS



Length	Width	Thickness	REF No non-sterile	REF No sterile
35 mm	13 mm	0.40 mm	281.1335.040	281.1335.040S



Length	Width	Thickness	REF No non-sterile	REF No sterile
39 mm	13 mm	0.40 mm	281.1339.040	281.1339.040S



Length	Width	Thickness	REF No non-sterile	REF No sterile
40 mm	15 mm	0.50 mm	281.1540.050	281.1540.050S



Length	Width	Thickness	REF No non-sterile	REF No sterile
41 mm	19 mm	0.60 mm	281.1941.060	281.1941.060S

**Saw Blade fits to Smith & Nephew DYONICS™ Saw REF 7205786 and 7205791**



Length	Width	Thickness	REF No non-sterile	REF No sterile
25 mm	6 mm	0.40 mm	282.0625.040	282.0625.040S



Length	Width	Thickness	REF No non-sterile	REF No sterile
25 mm	9 mm	0.40 mm	282.0925.040	282.0925.040S

**Saw Blade fits to Linvatec Hall®, MicroChoice®, Micro100™, Micro E®**



Length	Width	Thickness	REF No non-sterile	REF No sterile
18 mm	6 mm	0.40 mm	283.0618.040	283.0618.040S



Length	Width	Thickness	REF No non-sterile	REF No sterile
25 mm	6 mm	0.40 mm	283.0625.040	283.0625.040S
25 mm	6 mm	0.60 mm	283.0625.060	283.0625.060S



Length	Width	Thickness	REF No non-sterile	REF No sterile
25 mm	9 mm	0.40 mm	283.0925.040	283.0925.040S
25 mm	9 mm	0.60 mm	283.0925.060	283.0925.060S



Length	Width	Thickness	REF No non-sterile	REF No sterile
25 mm	9 mm	0.40 mm	283.0925.040F	283.0925.040FS



Length	Width	Thickness	REF No non-sterile	REF No sterile
35 mm	9 mm	0.60 mm	283.0935.060	283.0935.060S



Length	Width	Thickness	REF No non-sterile	REF No sterile
41 mm	9 mm	0.60 mm	283.0941.060	283.0941.060S



Length	Width	Thickness	REF No non-sterile	REF No sterile
25 mm	14 mm	0.60 mm	283.1425.060	283.1425.060S



Length	Width	Thickness	REF No non-sterile	REF No sterile
41 mm	14 mm	0.60 mm	283.1441.060	283.1441.060S



Length	Width	Thickness	REF No non-sterile	REF No sterile
41 mm	19 mm	0.60 mm	283.1941.060	283.1941.060S

**Saw Blade fits to Linvatec Hall®, MicroChoice®, Micro100™, Micro E®**



Length	Width	Thickness	REF No non-sterile	REF No sterile
14 mm	5 mm	0.40 mm	284.0514.040	284.0514.040S



Length	Width	Thickness	REF No non-sterile	REF No sterile
18 mm	5 mm	0.40 mm	284.0518.040	284.0518.040S



Length	Width	Thickness	REF No non-sterile	REF No sterile
25 mm	5 mm	0.40 mm	284.0525.040	284.0525.040S



Length	Width	Thickness	REF No non-sterile	REF No sterile
18 mm	9 mm	0.40 mm	284.0918.040	284.0918.040S



Length	Width	Thickness	REF No non-sterile	REF No sterile
35 mm	11 mm	0.40 mm	284.1135.040	284.1135.040S



## Saw Blade fits to Colibri (Synthes®)



Length	Width	Thickness	REF No non-sterile	REF No sterile
18 mm	5 mm	0.40 mm	285.0518.040F	285.0518.040FS



Length	Width	Thickness	REF No non-sterile	REF No sterile
25 mm	6 mm	0.60 mm	285.0625.060	285.0625.060S




Length	Width	Thickness	REF No non-sterile	REF No sterile
15 mm	7 mm	0.40 mm	285.0715.040F	285.0715.040FS



Length	Width	Thickness	REF No non-sterile	REF No sterile
31 mm	7 mm	0.40 mm	285.0731.040F	285.0731.040FS



Length	Width	Thickness	REF No non-sterile	REF No sterile
21 mm	9 mm	0.40 mm	285.0921.040F	285.0921.040FS




Length	Width	Thickness	REF No non-sterile	REF No sterile
22 mm	9 mm	0.40 mm	285.0922.040F	285.0922.040FS



Length	Width	Thickness	REF No non-sterile	REF No sterile
25 mm	10 mm	0.60 mm	285.1025.060	285.1025.060S



Length	Width	Thickness	REF No non-sterile	REF No sterile
50 mm	10 mm	0.60 mm	285.1050.060	285.1050.060S



Length	Width	Thickness	REF No non-sterile	REF No sterile
15 mm	11 mm	0.40 mm	285.1115.040F	285.1115.040FS



Length	Width	Thickness	REF No non-sterile	REF No sterile
27 mm	11 mm	0.40 mm	285.1127.040F	285.1127.040FS

## Saw Blade fits to Colibri (Synthes®)



Length	Width	Thickness	REF No non-sterile	REF No sterile
31 mm	11 mm	0.40 mm	285.1131.040F	285.1131.040FS



Length	Width	Thickness	REF No non-sterile	REF No sterile
40 mm	11 mm	0.50 mm	285.1140.050F	285.1140.050FS




Length	Width	Thickness	REF No non-sterile	REF No sterile
21 mm	13 mm	0.40 mm	285.1321.040F	285.1321.040FS



Length	Width	Thickness	REF No non-sterile	REF No sterile
25 mm	14 mm	0.60 mm	285.1425.060	285.1425.060S



Length	Width	Thickness	REF No non-sterile	REF No sterile
50 mm	14 mm	0.60 mm	285.1450.060	285.1450.060S
50 mm	14 mm	1.00 mm	285.1450.100	285.1450.100S



Length	Width	Thickness	REF No non-sterile	REF No sterile
15 mm	15 mm	0.40 mm	285.1515.040F	285.1515.040FS



Length	Width	Thickness	REF No non-sterile	REF No sterile
40 mm	15 mm	0.50 mm	285.1540.050F	285.1540.050FS



Length	Width	Thickness	REF No non-sterile	REF No sterile
50 mm	19 mm	0.60 mm	285.1950.060	285.1950.060S



Length	Width	Thickness	REF No non-sterile	REF No sterile
70 mm	19 mm	0.60 mm	285.1970.060	285.1970.060S

## Saw Blade fits to Colibri (Synthes®)



Length	Width	Thickness	REF No non-sterile	REF No sterile
50 mm	27 mm	0.60 mm	285.2750.060	285.2750.060S



Length	Width	Thickness	REF No non-sterile	REF No sterile
19 mm	10 mm	0.60 mm	285.1002.060	not available

## Saw Blade fits to EPEN (Synthes®)



Length	Width	Thickness	REF No non-sterile	REF No sterile
18 mm	3 mm	0.40 mm	286.0318.040F	286.0318.040FS



Length	Width	Thickness	REF No non-sterile	REF No sterile
18 mm	5 mm	0.40 mm	286.0518.040F	286.0518.040FS



Length	Width	Thickness	REF No non-sterile	REF No sterile
27 mm	6 mm	0.60 mm	286.0627.060	286.0627.060S



Length	Width	Thickness	REF No non-sterile	REF No sterile
15 mm	7 mm	0.40 mm	286.0715.040F	286.0715.040FS



Length	Width	Thickness	REF No non-sterile	REF No sterile
31 mm	7 mm	0.40 mm	286.0731.040F	286.0731.040FS



Length	Width	Thickness	REF No non-sterile	REF No sterile
22 mm	9 mm	0.40 mm	286.0922.040F	286.0922.040FS



Length	Width	Thickness	REF No non-sterile	REF No sterile
27 mm	10 mm	0.60 mm	286.1027.060	286.1027.060S



Length	Width	Thickness	REF No non-sterile	REF No sterile
15 mm	11 mm	0.40 mm	286.1115.040F	286.1115.040FS



Length	Width	Thickness	REF No non-sterile	REF No sterile
31 mm	11 mm	0.40 mm	286.1131.040F	286.1131.040FS



Length	Width	Thickness	REF No non-sterile	REF No sterile
22 mm	13 mm	0.40 mm	286.1322.040F	286.1322.040FS

## Saw Blade fits to EPEN (Synthes®)



Length	Width	Thickness	REF No non-sterile	REF No sterile
27 mm	14 mm	<b>0.60 mm</b>	<b>286.1427.060</b>	<b>286.1427.060S</b>



Length	Width	Thickness	REF No non-sterile	REF No sterile
15 mm	15 mm	<b>0.40 mm</b>	<b>286.1515.040F</b>	<b>286.1515.040FS</b>



Length	Width	Thickness	REF No non-sterile	REF No sterile
41 mm	19 mm	<b>0.60 mm</b>	<b>286.1941.060</b>	<b>286.1941.060S</b>

## Saw Blade fits to Arthrex



Length	Width	Thickness	REF No non-sterile	REF No sterile
18 mm	5 mm	0.40 mm	287.0518.040F	287.0518.040FS



Length	Width	Thickness	REF No non-sterile	REF No sterile
25 mm	6 mm	0.60 mm	287.0625.060F	287.0625.060FS



Length	Width	Thickness	REF No non-sterile	REF No sterile
25 mm	9 mm	0.40 mm	287.0925.040	287.0925.040S



Length	Width	Thickness	REF No non-sterile	REF No sterile
25 mm	9 mm	0.40 mm	287.0925.040F	287.0925.040FS



Length	Width	Thickness	REF No non-sterile	REF No sterile
35 mm	9 mm	0.60 mm	287.0935.060	287.0935.060S



Length	Width	Thickness	REF No non-sterile	REF No sterile
41 mm	9 mm	0.60 mm	287.0941.060	287.0941.060S




Length	Width	Thickness	REF No non-sterile	REF No sterile
41 mm	9 mm	0.60 mm	287.0941.060F	287.0941.060FS




Length	Width	Thickness	REF No non-sterile	REF No sterile
41 mm	14 mm	0.60 mm	287.1441.060	287.1441.060S


## Saw Blade fits to Aesculap® Acculan 3Ti




Length	Width	Thickness	REF No non-sterile	REF No sterile
25 mm	6 mm	0.60 mm	289.0625.060	289.0625.060S




Length	Width	Thickness	REF No non-sterile	REF No sterile
31 mm	7 mm	0.40 mm	289.0731.040F	289.0731.040FS



Length	Width	Thickness	REF No non-sterile	REF No sterile
21 mm	9 mm	0.40 mm	289.0921.040F	289.0921.040FS




Length	Width	Thickness	REF No non-sterile	REF No sterile
25 mm	9 mm	0.40 mm	289.0925.040	289.0925.040S



Length	Width	Thickness	REF No non-sterile	REF No sterile
25 mm	10 mm	0.60 mm	289.1025.060	289.1025.060S



Length	Width	Thickness	REF No non-sterile	REF No sterile
50 mm	10 mm	0.60 mm	289.1050.060	289.1050.060S



Length	Width	Thickness	REF No non-sterile	REF No sterile
25 mm	14 mm	0.60 mm	289.1425.060	289.1425.060S

## Saw Blade fits to Aesculap® Acculan 3Ti



Length	Width	Thickness	REF No non-sterile	REF No sterile
40 mm	15 mm	0.50 mm	289.1540.050	289.1540.050S



Length	Width	Thickness	REF No non-sterile	REF No sterile
50 mm	19 mm	0.60 mm	289.1950.060	289.1950.060S