

GREETING

We nelis, are not just satisfied as a certified manufacturer, steadily studying and developing and developing our technique under the company's motto "creation" to provide all surgeons the most convenient and the safest medical instruments for their patients. All staff will spare no efforts to make our desires come true!

We have been devoting its efforts to broaden experiences and research in bio medical engineering to become a leading stent manufacturer.

All members in nelis, with ceaseless challenge and creative mind, will do our best to satisfy the customers needs and also maintain our effort to contribute to healthier life of people.

Thanks for your support and encouragement.

CEO *Onam Seo*



CERTIFICATIONS

	APPROVED ITEM	LICENCE NO.	AUTHORITY	APPROVAL DATE
GMP	RETRACTOR, TROCAR, BIOPSY PROCEDURE KIT	1483	KFDA	13.04.2010
ISO 9001:2008	SINGLE PORT SPECIMEN BAG	KR10 / 1466	SGS	21.06.2010
ISO 13485:2003		KR10 / 1467		
CE 0120		KR10 / 01468		
FDA	GLOVE PORT	K141715	FOOD AND DRUG ADMINISTRATION	25.03.2015
	ENDO KEEPER	K152598		17.03.2016

Company's information

Address	#1005, 201-dong, Bucheon Techno Park Ssangyong 3-cha, Samjeong-dong, Bucheon-si, Gyeonggi-do, 14449, Korea
Tel	+82 32 624 1697
Fax	+82 32 624 1699
Email	nelismed@daum.net

GLOVE PORT

This is the real single port.



GLOVE PORT A version

Glove port A version can be open during operation.
It's easy to extract with surgical specimen during operation

Glove port A Simple characteristics

- 1) Able to open during operation
- 2) Robotic surgery
- 3) Separation of wound retractor and port
- 4) Small scar
- 5) Easy insert ring design
- 6) 3mm port ~ 15mm port
- 7) Separation of gas inlet and outlet

Applicable operation

Operation	Model no.
Appendectomy	440, 431, 431-AS, 431-AM
Gall Bladder	431, 431-AS, 440, 440-AS
Uterine Myoma	330, 431, 431-AS, 330-AS
Hysterectomy	422, 422-AS, 422-AM
Hysteromyomectomy with morcellator	P422-15, 422-AS-15, 422-AM-15
Cholecystectomy	440, 431, 431-AS, 431-AM
Ovariectomy	330, 431, 330-AS, 431-AS
Urinary Tract Stone Removal	P431, 431-AS, 431-AM
Endometriosis	P312, 321-AS, 330-AS
Paediatric surgery	431-AT, 440-AT, 4220-AT
Urology	431, 422, 431-AS, 431-AM,
Colorectal surgery	431-AS, 422-AS, 431-AM, 422-AM,
Gastrectomy	422-AM, 431-AM

*OTHERS : Robot surgery, Ectopic pregnancy, Fibroids, Hepatectomy and so on



ENDO KEEPER CG series

Endo Keeper CG series enables surgical instruments through the port (5-12mm) and hands to be put in the body and take out large specimens.

Endo Keeper CG version simple characteristics

- Easy insert with two step design
- Adjust the height to the thicknesses of the patients' abdominal
- High elasticity
- Easy removal ring
- Rolling the sheet helps to open incision easily
- Can be opened during operation
- Size of ports 5mm~12mm
- Separation gas inlet
- CO₂ injection is possible without the help of other instrument

Endo Keeper Model List & Measurement

	Length	Diameter	Incision
X-SAMLL	130mm	41mm	20mm
SMALL-S	185mm	60mm	20-40mm
SMALL-NS			
MEDIUM-S	185mm	95mm	25-60mm
MEDIUM-NS			
LARGE-N	190mm	125mm	50-90mm
XX-LARGE	210mm	195mm	90-150mm
SMALL-NS-C	265mm	60mm	20-40mm
SMALL-NS-CG			
SMALL-S-C	265mm	60mm	20-40mm
SMALL-S-CG			
MEDIUM-NS-C	265mm	95mm	25-60mm
MEDIUM-NS-CG			
MEDIUM-S-C	265mm	95mm	25-60mm
MEDIUM-S-CG			



**SINGLE PORT
LAPAROSCOPIC SURGERY**

WOUND RETRACTOR