



Smart Scale - Uroflowmetry

NON INVASIVE DIAGNOSTICS.
SEVERAL SOLUTIONS FOR DIFFERENT NEEDS.



Uroflowmetry B Series

URODYN+ FLOWMETER FEATURES

- Straightforward for fast & easy uroflow reading.
- Wireless system maximizes patient privacy to make results more natural.
- Real-time flow display on-screen.
- Intuitive touchscreen software.
- Quick guides for setup & training.
- Rapid report printout.



SPINNING DISK

- Unique technology for the most accurate & reliable uroflow available.
- Potent software filter to screen out artifacts.

TECHNICAL SPECIFICATIONS

Items	Details	Items	Details
Range	Urine flow rate: 0~100 mL/s, error $\pm 1\%$ or 1 mL/s, whichever is greater	Urine volume	0~1100 ml, error $\pm 1\%$ or 1ml, whichever is greater
Size	Urine collector: 270*290*350mm	Urination time	0s~240s; measurement error: $\pm 1s$
Weight	Urine collector: 3.0KG	Maximum volume	1100ml

TECHNICAL SPECIFICATIONS

- 8" Touchscreen Size.
- Integrated Thermal Printer for Report Printing.
- Optional Rechargeable Battery.
- Optional Wheeled Floor Stand.

