

SUCTION CONNECTING TUBE CONNECTING TUBE



Suction Connecting Tube™ GS-5017

- Suction connecting tube with Molded Connectors.
- Flexible connectors ensure easy and secure connection.
- Sterile, Individually packed.
- Box of 10/ Master Box of 100.

Suction Connecting Tube™ Integral Funnel Ends GS-5019

- Suction connecting tube with Integral Funnel ends.
- Easy connection to suction cannisters on one end, and Yankaur suction handle/ tip at the other end.
- The integral funnel connectors can be modified by cutting with scissors for extension, if required.
- Sterile, Individually packed.
- Box of 10/Master Box of 100.
- Non-Sterile rolls available with bubble at each 100 cm.

SPECIFICATIONS	
Tube Length	I.D/ O.D
GS-5017	
200, 250, 300 cm	Tube: 6.0 / 9.0 mm
GS-5019	
180, 220, 270, 300 cm	Tube: 7.0 / 10 mm Ends: 10 / 15 mm
Optional: 30 mtr roll (Non-Sterile)	Tube: 7.0 / 10 mm Ends: 10 / 15 mm

CONDUCTIVE SUCTION TUBE CONNECTING TUBE



Conductive Suction Tube™ GS-5101

- Suction connecting tube with conductive line to prevent static charge.
- Provides easy connection to suction cannisters on one end and yankaur suction handle/ tip at the other end.
- The integral funnel connectors can be modified by cutting with scissors for fixation/ connection as per requirement.
- Clear, non-coiling tube's wall thickness prevents it from collapsing when under high negative pressure.
- Sterile, Individually packed in a peelable pouch.
- Box of 10/ Master Box of 100 pieces.
- Non-Sterile rolls available with bubble at each 100 cm.

SPECIFICATIONS	
Tube Length	I.D/ O.D
180, 220, 270, 300 cm	Tube: 5.0 / 7.0 mm Ends: 10 / 13.5 mm
Optional: 30 mtr roll (Non-Sterile)	Tube: 5.0 / 7.0 mm Ends: 10 / 13.5 mm