

**CMOS+LED**

**FLEXIBLE VIDEO INTUBATION LARYNGOSCOPE**





## ● Multi-function Buttons

Capture, record, freeze, zoom in/out.



\* For QG-3032, QG-3048 and QG-3052 only.

## ● Detachable Battery

Always ready for use.



## ● Portable Type

Compatible with portable monitor CV-1080/HV-3101.



## ● Handheld Type

Compatible with CV-1030.




## CMOS + LED Technology

CMOS image sensor delivers a high quality image.  
No fiber: durable and easy maintenance.



## Application

Gold standard for difficult airway management:  
intubation, positioning, maintenance and  
interventional airway treatment.



Model	Insertion Tube(mm)	Working Channel(mm)	Working Length(mm)	Deflection Angle	View Angle	Depth of View	Image Processor
QG-2128	2.8	-	600	180°/130°	90°	3-50mm	CV-1030 CV-1080 HV-3101
QG-2138	3.8	-			90°		
QG-2152	5.2	2.2			90°		
QG-3032	3.1	1.2			120°		
QG-3048	3.8	1.5			120°		
QG-3052	5.2	2.6			90°		