



◆ Ligating Clips


- Used for human body blood vessels and coelom organization ligation and closing during operation.



LC-1-6



LC-2-6




LC-3-6

- ▶ Safe polymer materials, good biological compatibility.
- ▶ No affecting on X-ray, CT, MRI and other imaging studies.
- ▶ Safety lock, arc, elasticity, and prevent slippery designs. Rapid ligation in operation, safe and reliable results.

Recommendation Specification

Code	Size	Length of Clips (mm)	Scope of Closure Diameter (mm)
LC-1-6	ML	10	3~10
LC-2-6	L	13	5~13
LC-3-6	XL	17	7~16

◆ Applier



Open Applier

Code	Length	Angle	color	Matched clip size
CA C20*9S	20CM	9°	Green	ML
CA C27*9S	27CM	9°	Green	ML
CA C20*70S	20CM	70°	Green	ML
CA C27*70S	27CM	70°	Green	ML
CA C20*9M	20CM	9°	Purple	L
CA C27*9M	27CM	9°	Purple	L
CA C20*70M	20CM	70°	Purple	L
CA C27*70M	27CM	70°	Purple	L
CA C20*9L	20CM	9°	Gold	XL
CA C27*9L	27CM	9°	Gold	XL
CA C20*70L	20CM	70°	Gold	XL
CA C27*70L	27CM	70°	Gold	XL



Laparoscopic Applier

Code	Size	Color	Matched clip size
CA B 5-S	Φ5mm x 330mm	Green	ML
CA B 10-M	Φ10mm x 330mm	Purple	L
CA B 10-L	Φ10mm x 330mm	Golden	XL



Remover

Code	Size	Matched Clip Size
CR E-SM	Φ5mm x 330mm	ML, L
CR E-ML	Φ10mm x 330mm	L, XL