

TROCAR FILTERING SYSTEM



TROPIAN

Good grip : Ergonomic designed trocar palm.

Stable position hold : Ribbed Cannula can be firmly fixed on the surgical part on the abdominal wall.

Dual Seal shields : Minimum leakage CO2 gas during operation.

Bladless trocar tip : Abdominal wall rapid reconstruction.

Clean site : Transparent Cannula enables visual control of every surgical step.

Responsible for the health of surgeons and assistants : Equipped with filter for smoke evacuation.
(2 model : ZT-12F, ZT-22F)

Integrated converter : Tropan (ZT-30, ZT-40) is possible to use 5-12mm diameter instrument without converter.



Dual Seal shields

Minimum leakage CO2 gas during operation



Responsible for the health of surgeons and assistants

Equipped with filter for smoke On / Off possible



Bladless trocar tip

Abodminal wall rapid reconstruction



Transparent Cannula enables visual control of every surgical step.





TROPIAN TECH Co., Ltd.

Tropian Specification Sheet

Tropian Trocar System

Product	Code	Description	Length	Quantity
	ZT-10	Bladeless Trocar, Ø5mm	70mm	10/box
	ZT-20	Bladeless Trocar, Ø5mm	100mm	10/box
	ZT-30	Bladeless Trocar, Ø11mm	100mm	10/box
	ZT-40	Bladeless Trocar, Ø12mm	100mm	10/box
	ZTC-10	Thoracoscopic Bladeless Trocar, Ø5mm	75mm	10/box
	ZTC-30	Thoracoscopic Bladeless Trocar, Ø11mm	75mm	10/box
	ZTC-40	Thoracoscopic Bladeless Trocar, Ø12mm	75mm	10/box

Tropian Trocar Filtering System

Product	Code	Description	Length	Quantity
	ZT-12F	Bladeless Trocar with filter, Ø5mm	70mm	10/box
	ZT-22F	Bladeless Trocar with filter, Ø5mm	100mm	10/box



TROPIAN TECH Co., Ltd.

#302 3F ISU Building

9, song I-ro 14gil, Songpa-ku, Seoul, Korea

Phone +82 2 402 6315 **Fax** +82 2 402 6310

www.tropiantech.com

office@tropiantech.com