

# EZLIG<sup>®</sup>

## MULTI - BAND LIGATOR

Clearer field & Easier to operate



CE 0123

<http://www.bjzksk.com>

# EZLIG<sup>®</sup>

## MULTI - BAND LIGATOR

### INDICATIONS


Used for ligating esophageal varices combined with Endoscope.

### TECHNICAL AND FUNCTIONAL FEATURES


- Post installation of the bands makes the field clear
- Precise ligation, obvious release resistance, easy to operate
- Multi-bands design achieves to treat several points in one operation
- With good biocompatibility and physical property
- Humanized handle design make operation more easy and comfortable
- Bands release reminds by feeling
- Different code color rings for quick identification of the last shot
- Square shape cross section design make ligation stronger and safer



### Excellent filed & aspiration effect



### Single & multiple variceal ligation



### SPECIFICATIONS

Model	Band No.	Trigger Cord
BL18-7A	7	0.2mm
BL18-9A	9	0.2mm

Beijing ZKSK Technology Co., Ltd

Building 9, 6 & No.6 Yuan Hengye North 7th street,  
Yongle Economic Development Zone, Tongzhou District,  
Beijing 101105, China

Tel: 86 10 57272222

Fax: 86 10 63777521

E-mail: bjzksk3@163.com

<http://www.bjzksk.com>

EC REP Prolinx GmbH

Brehmstr. 56, 40239 Duesseldorf, Germany

CE 0123