

RTA-19 Incubator																						
	Features I. Compact, dry heater. II. Easy use with Chemistry Analyzer III. Incubates tubes at 25 °C/30 °C/37 °C/45 °C IV. LED temperature display V. Two kinds of tube, up to 24 tubes incubated																					
	System Specifications <table> <tr> <td>Heat up mode</td><td>Metalline heater</td></tr> <tr> <td>Temperature option</td><td>25 °C/30 °C/37 °C/45 °C</td></tr> <tr> <td>Temperature precision</td><td>±0.5 °C</td></tr> <tr> <td>Temperature Difference between holes</td><td>±0.3 °C</td></tr> <tr> <td>Heat up speed</td><td>25 °C~37 °C within 6 minutes (room temperature 25 °C)</td></tr> <tr> <td>Temperature Display</td><td>Digital LED</td></tr> <tr> <td>Hole specification</td><td>Φ13mmX20 、 Φ15mmX4 ; hole depth: 33mm</td></tr> <tr> <td>Heat up power</td><td>70W</td></tr> <tr> <td>Dimensions L × W × H(mm)</td><td>272mm x 203mm x 188mm</td></tr> <tr> <td>Net Weight</td><td>2.4kg</td></tr> <tr> <td>Power supply</td><td>AC 85 ~ 265V 50 ~ 60HZ</td></tr> </table>	Heat up mode	Metalline heater	Temperature option	25 °C/30 °C/37 °C/45 °C	Temperature precision	±0.5 °C	Temperature Difference between holes	±0.3 °C	Heat up speed	25 °C~37 °C within 6 minutes (room temperature 25 °C)	Temperature Display	Digital LED	Hole specification	Φ13mmX20 、 Φ15mmX4 ; hole depth: 33mm	Heat up power	70W	Dimensions L × W × H(mm)	272mm x 203mm x 188mm	Net Weight	2.4kg	Power supply
Heat up mode	Metalline heater																					
Temperature option	25 °C/30 °C/37 °C/45 °C																					
Temperature precision	±0.5 °C																					
Temperature Difference between holes	±0.3 °C																					
Heat up speed	25 °C~37 °C within 6 minutes (room temperature 25 °C)																					
Temperature Display	Digital LED																					
Hole specification	Φ13mmX20 、 Φ15mmX4 ; hole depth: 33mm																					
Heat up power	70W																					
Dimensions L × W × H(mm)	272mm x 203mm x 188mm																					
Net Weight	2.4kg																					
Power supply	AC 85 ~ 265V 50 ~ 60HZ																					



Rayto Life and Analytical Sciences Co.,Ltd.
 A/6F,Technology Building,Industry Road,Shekou,Shenzhen P.R.China
 Tel:+86 755 86168182 Fax:+86 755 86168796
 Website: www.rayto.com