

UP-897MD

Black & White
A6 Video Graphic Printer

The Sony UP-897MD black and white video graphic printer is designed specifically for use with medical diagnostic equipment, such as ultrasound systems. The UP-897MD packs power in a small chassis. High quality, photo-like prints are output in just approximately 2.0 seconds. The UP-897MD is extremely small and light weight. For additional convenience the UP-897MD adds a jog dial to the front panel for easy operation.



UP-897MD

- Ideal for ultrasound printing
- Smaller and lighter by 12% compared to the 895 series printer
- User-friendly controls on front panel with addition of jog dial
- High quality, photo-like prints in approximately 2 seconds
- Uses the same media as the UP-895 series printers: UPP-110S/5, UPP-110HD, UPP-110HG

SONY®



Sony Medical

UP-897MD

Specifications

UP-897MD

Printing Method	Direct Thermal Printing
Thermal Head	Thin film thermal 1280 dots
Resolution	325 DPI
Gradations	True 256 gray balance
Effective Pixels	EIA 1280 dot x 508 lines CCIR: 1280 dot x 608 lines
Memory	10 Frames
Print Speed	Approximately 2.0 seconds (High Speed Mode)
Input/Output	Video In: BNC, CCIR or EIA composite video signal, 1.0 Vp-p, 75Ω / high impedance selectable Video Out: BNC, Loop-through or E to E selectable
Remote Terminal	Stereo Mini-jack
Power Requirements	100 to 240 V AC, 50/60Hz, 1.5A to 0.8A
Dimensions	154mm x 88mm x 240mm (WHD) (6 1/8 x 3 1/2 x 9 1/2 inches – WHD)
Weight	Approximately 2.6kg (5lb. 12 oz)
Safety Standards	UL60950-1 Listed + CSA C22.2 No. 60950-1 + UL60601-1 Listed + CSA No. 60601-1
Supplied Accessories	AC Power Cord Operation Manual (CD ROM) BNC Video Connecting Cable Thermal Head Cleaning Sheet UPP-110HG Media (1Roll)
Media (Not Included)	UPP-110HG High Gloss Paper UPP-110HD High Density Paper UPP-110S/5 Standard Paper
Optional Accessories	FS-20 Wired foot switch RM-91 Remote thumb button



UP-897MD Back



Print Media: UPP-110S/5 (Standard media),
UPP-110HD (High density media),
UPP-110HG (Premium high-gloss media)

SONY

Sony Electronics Inc.
1 Sony Drive, Park Ridge, NJ 07656

www.sony.com/medical