

UA-66

Urine Analyzer



- Cold light source with high luminosity, long shelf life
- Innovative light sensor, high sensitivity
- Automatic calibration improves accuracy
- Flag abnormal results
- Standard RS232 interface for connection with external computer
- Storage of 1000 testing results
- Multi-formats units
- Built in thermal printer
- Large LCD Display

UA-66

Urine Analyzer

Specifications

Test items	Urobilinogen, Bilirubin, Ketone, Blood, Protein, Nitrite, Leukocytes, Glucose, Specific Gravity, Ascorbic Acid, pH
Wavelength of monochromatic light	557nm, 620nm, 656nm
Principle	Reflectance photometer
Strips	11 items
Throughput	120 tests/hour or 60 tests/hour
Display	Large LCD display, Resolution: 240x128
Memory	1000 results
Printer	Built in thermal printer, paper width: 57mm
Interface	Standard RS232
Working environment	Temperature 5°C~40°C Humidity 30%~85%
Power supply	220V~, 50Hz 110V~, 50Hz
Power	40VA
Dimension	100mm(H)x360mm(W)x300mm(D)
Weight	3.8kg



U-11

Urine Strips

Urobilinogen, Bilirubin, Ketone, Blood, Protein, Nitrite, Leukocytes, Glucose, Specific Gravity, Ascorbic Acid, pH

Accurate Results

High sensitivity and specificity

Immediate Diagnosis

Result in one minute

Convenient

Simple procedure

