

ENDO



ENDO EI.EC.v7

EI.LS.v7

Video Endoscope System & Light Source

<https://endo.id>

Innovation rewrites the standard. With that in mind, ENDO Indonesia presents you the ENDO EI.EC.v7, a brand new high definition endoscopy system designed specifically for screening and diagnosis with an excellent budget friendliness.





CMOS

With the CMOS image sensor, ENDO EI.EC.v7 gives high quality images with low noise levels, which leads to improved clarity and optimized vision in not only routine procedures but also in complex interventions.



Smart Workflow

To increase efficiency, patient data management is accessible with the press of button. The smart workflow of the built-in workstation of ENDO EI.EC.v7 helps doctors to edit reports, record videos and manage patient files with ease.



Clean and Clear

The waterjet function of the scopes helps to maintain a clear endoscopic view by washing away blood or mucous.



ENDO EI.EC.v7 Video Endoscope System

HD video outputs	DVI,VGA
Standard video outputs	S-Video, CVBS
HD video recording	Built-in workstation, timely record and review
Zoom / PBP	Yes / Yes
Dimension (L x W x H / weight)	455mm x 370mm x 124mm / 9.5 kg

ENDO EI.LS.v7 Light Source


Main lamp	50W LED
Transparency function	7 seconds maximum brightness flicker
Lifetime indicator	Yes
Hot plugging	Endoscope switch without powering off
Dimension (L x W x H / weight)	460mm x 370mm x 154mm / 11kg

430 Series Scopes

Parameter	Gastroscope		Colonoscope			
	EG-430	EG-430L	EC-430	EC-430T	EC-430L	EC-430L/T
Depth of view (mm)	3-100		3-100			
Field of view	140°		140°			
Distal end diameter (mm)	9.3	9.8	12	12	12.9	12.9
Bending section diameter(mm)	9.3	9.8	12.5	12.5	12.9	12.9
Instrument channel diameter (mm)	2.8	3.2	3.8		4.2	
Angulations	UP 210°/DOWN 90° LEFT 100°/RIGHT 100°		UP 180°/DOWN 180°		LEFT 160°/RIGHT 160°	
Working length (mm)	1050		1350	1700	1350	1700
Total length (mm)	1350		1700	1950	1700	1950


ENDO

 Jalan Raya Menganti 14 Kedurus, Surabaya 60223
Jawa Timur, Indonesia

 62-31-7673636

 62-31-7673737

 0800 . 177 . ENDO (3636)

 info@endo.id

 <https://endo.id>