

# ENDO



## **ENDO ANM3** **Anesthesia Machine**

<https://endo.id>



## Featured

- CE marked, meet EU clinical requirements
- Contemporary and ergonomic designed for ease of use
- Pull handle
- DIN rail at both side of machine allows for flexibility and ease of mounting patient monitor and other expansion medical devices
- Optional Integrated auxiliary O<sub>2</sub> flowmeter and suction control
- AGSS (Anesthesia Gas Scavenging System) is optional.
- Electronic displayed flowmeter
- Applicable to Infant, pediatric and adult.
- Driven gas can be Oxygen or Air. save running cost of Oxygen.

## Integrated CO<sub>2</sub> Absorber System

- Highly integrated design - fewer parts and connections reduces potential for leaks and mismatches.
- Ease of disassembly (no tools required).
- Ease for cleaning, autoclavable in 134°C to avoid cross-infection (except for O<sub>2</sub> sensor and gauge), latex free.
- One step bag/vent, switch turns the ventilator on/off.
- Automatic CO<sub>2</sub> absorber by-pass.
- ACGO.

## Ventilator

- 12 inch touch screen for easy review.
- Touchscreen + navigation wheel+ key pad for simple operation.
- "Choose-set-confirm" to realize function and avoid mis-operation.
- One key presses to total Standby mode: end case
- Easy of use, full range ventilator with mechanical ventilation modes, spontaneous ventilation modes and electronic, support a variety of patients
- Powerful extended function: Optional plug-in anesthesia agent monitoring CO<sub>2</sub> module, and SpO<sub>2</sub> module. With parameters setting/display and graphic display in the screen.

- World leading brands of flow sensors, differential pressure sensors, pressure regulating valves and proportional valves to ensure accurate ventilation.
- Tidal volume compensation provides accurate delivery down to 20 ml
- Smart alarm notification direct user to specific problems and affected parameters
- Pressure and volume waveforms on a breath-by-breath basis
- Audio and visual alarms with three levels.

## Vaporizers

### Anesthetic vapour delivery modules

- The highest quality and reliable vaporizer.
- Available for Halothane, Isoflurane, Enflurane and Sevoflurane.
- Service-free: Calibration every 10 years (5 years for Halothane).
- Superb performance with automatic compensation for temperature, pressure and flow.
- Selectatec compatible interlock system.