

ENDO

ENDO Baby Incubator EI.BI2



Spesifikasi

Power supply	AC110/220V 60/50Hz	Power input	650 VA
Kondisi Pengoperasian			
Suhu	20~30°	Keseragaman suhu kasur	≤0.8°C
Kelembaban relative	30~75%	Tingkat kebisingan intern	<55dB (A)
Kecepatan aliran udara	<0.3 m/detik	Sudut kemiringan	±12°
Rentang kontrol suhu udara	25°C ~ 37°C 37.1°C ~ 38°C	Ukuran kasur	65cm (L) x 37cm (W)
Kisaran kendali suhu kulit	32°C ~ 37°C	Kapasitas Palung	1000 mL
Akurasi sensor suhu kulit	±0.2°C	Penyaring udara	0.5 µm
Fluktuasi suhu	±0.5°C	Dimensi Setiap unit dikemas dalam satu kotak polywood (unit kontrol dan alas)	126 x 71 x 85 cm
Berat	Net 80Kg Gross 100Kg		