

ENDO



ENDO DRS

Digital Radiography System

<https://endo.id>

ENDO DRS

Digital Radiography System

Anywhere You Can Get

Incredibly easy to move, and simple to operate, ENDO DRS is ideal for various radiography applications. High flexibility through the use combined with high-quality images at 5.6 kW output allows you to get quick and dependable X-ray imaging anywhere you need it.



5.6 kW high power



Internal Li-ion battery



Dual touch screen control



Rotatable X-ray unit



Storage bin with locking system



Ergonomic handle



Lightweight and small occupation



Wireless detector

Smooth Operation, Enhanced Efficiency

Imaging without limits

- ENDO DRS supports working both from mains and internal battery power supply.
- High performance of the lithium-ion battery in long standby time and short charging time guarantee a stable power supply.

Freely positioning

- Cable-free connection among X-ray unit, image workstation and wireless detector permits positioning without constraints.
- Rotatable X-ray unit keeps maneuverable even in restricted spaces, making it suitable for bedside imaging.
- Benefiting from smooth positioning, ENDO DRS can shorten patients' wait times and maximize your work efficiency.

All-in-one solution

- Integrated with an image workstation, ENDO DRS brings a fully digital workflow as required.
- Large touch screen enables smart control of the complete system as well as great visibility of details.

Simple to operate

- Equipped with multiple preset APRs, an 8-inch touch screen with user-friendly interface results in extremely convenient operation.


Transport with reliability


- Lightweight body combined with unobstructed and clear line of sight provides safe and easy movement.
- The large wheels ensure easy access to all terrain and won't tip over on uneven roads.
- Storage bin with a flip-up cover for detectors can greatly protect from dust and splash while in transport or storage.



Intelligent Workstation

- Fully digital workflow via DICOM 3.0 standard.
- Simple and intuitive user interface speeds up imaging workflow.
- Images can be sent directly to your PACS, and/or imager quickly, which maximizes the efficiency of the clinical diagnosis process.





Fast imaging



Frequency division



Noise reduction



Imaging enhancement





LUT curve



Balanced processing

Clinical images



ENDO

Jalan Raya Menganti 14 Kedurus, Surabaya 60223
Jawa Timur, Indonesia

62-31-7673636

62-31-7673737

0800 . 177 . ENDO (3636)

info@endo.id

<https://endo.id>