

# ENDO

## Camera System Integrated With Light Source ENDO EI.EC14



### Spesifikasi

#### Data Umum

Tegangan saluran	100-240 VAC
Persyaratan daya	95VA
Sekering listrik	Sekering, 5x20mm, inert 3.15AH, 250V dengan kapasitas switching yang tinggi (Ia = 1500A) Unit control: 295 x 128 x 278 mm
Dimensi (LxTxD)	Camera head tanpa tombol: ± 036.5 x 65 mm (tanpa grommet) Camera head tanpa tombol: ±038 x 79 mm (tanpa grommet)
Berat	Unit control: 5kg Camera head tanpa tombol: 150g Camera head dengan tombol: 180g

## Light Source

**Lightsource:**  
LED Satuan LED; sebanding 50W

## Kamera

Sensor gambar 1/3" MOS

System TV Dihubungkan pada 50/60 Hz

Resolusi TV 1920 x 1080 pixel (Full HD)

Resolusi gambar >900 TV lines; horizontal

Rasio indeks 16:09

Sensitivitas minimal <3lux (1/1.4)

## Koneksi

Output video digital  
Koneksi lainnya

- 1x HDMI/DVI
- 1x PS2 (keyboard)
- 1x Remote IN (footswitch)
- 1x Remote OUT (Videoprinter, Videorecorder, ...)
- 1x RJ45 (Service interface)

## Klasifikasi / Kesesuaian

Kelas keselamatan 1  
Bagian aplikasi BF