

ENDO

ENDO 3 in 1 GCU Meter EI.GCU



3 in 1, One test Three result

Glucose, Uric Acid, Cholesterol

Jangkauan deteksi luas

Strip

ENDO Multi Test Strip (EI.MTS)

Spesifikasi

Rentang pengukuran	GLU: 0.6~33.3mmol/L (10.8~600mg/dL) UA: 119~1190 μ mol/L (2~20mg/dL) TC: 2.59~12.93mmol/L (100~500mg/dL)
Tipe sampel	Darah kapiler dan vena utuh
Volume sample	35-40 μ L
Hasil waktu	\leq 120 detik
Penyimpanan	Sampai dengan 500 hasil sampel
Dimensi	142mm (L) x 70mm (W) x 26.2mm (H)
Berat	127g
Data output	Mini USB
Daya	DC 4.5V (3 AAA baterai)