

ENDO

ENDO Glukometer EI.GM.2




2 in 1, GLU + UA

User friendly, rentang volume sample kecil

Strip

ENDO Blood Glucose Test Strip (EI.GTS.2)

ENDO Uric Acid Test Strip (EI.UTS)



Three-electrode Strip

Spesifikasi

Rentang pengukuran	Uric Acid: 150-1200 μ mol/L Blood Glucose: 1.1-33.3mmol/L
Tipe sampel	Darah kapiler atau vena utuh
Volume sample	Uric Acid: 1.5 μ L Blood Glucose: 1.0 μ L
Hasil waktu	Uric Acid: 25 detik Blood Glucose: 5 detik
Penyimpanan	Sampai dengan 400 hasil sampel
Dimensi	104mm (L) x 60mm (W) x 21mm (H)
Berat	65g (dengan baterai)
Daya baterai	Setidaknya 1.000 tes
Daya	DC 3V (1 CR2032 baterai)