

ENDO

ENDO BLOOD GLUCOSE TEST STRIP EI.GTS



Cocok digunakan dengan
ENDO Glukometer (EI.GM)

SPESIFIKASI TEKNIS

Volume sample darah	≥ 0.8μL
Waktu pengujian	5 detik
Rentang pengukuran	0.5 - 33.5 mmol/L (9.0 – 600.0 mg/dL)
Suhu	23±4 °C (optimum)
Kelembapan	≤ 90%
Tipe sample	Darah kapiler
Kemasan	25 Strip/botol;
	320 botol/karton