

ENDO

ENDO Automated Hematology Analyzer EI.HA.10



- **5-Part-Diff Hematology Analyzer**
- **2D Grams**
Two 2D histograms: RBC and PLT
Two 2D scattergrams: 5-part differential
Eosinophi and neutrophil
- Layar sentuh sederhana dengan ikon bergambar yang ramah pengguna
- **Penyimpanan hingga 200.000**
- **Laser Scatter**
Multi-angle scatter separation (MASS)
Stabil dan long-lasting
- **Diluting unit**

Spesifikasi	
Display	10.4-inch LCD
Dimensi	530mm(L)×270mm(W)×412mm(H)
Berat	26.5Kg
Power supply	AC 100V~240V
Frekuensi	50/60Hz
Power	250VA
Parameter	34 parameter, 2 histograms, 2 scatter grams
Volume aspirasi	Whole blood 20 µL, pre-diluted 20µL
Input/output	Built-in thermal printer, support serial port, LAN port, USB port