

ENDO

ENDO Examination Headlamp EI.HL.2



Spesifikasi

Light source	Single LED
LED power	5W
Umur rentang LED	≥ 50.000 jam
Iluminasi	55000-65000 lux
Work distance	200-500mm
Light spot	Variable dari D=45mm hingga D=110mm
Suhu warna	3000-6500K
Power supply	AC 90-240V/DC5V
Baterai	Single (Standard), Dual (Optional)
Kapasitas baterai	1800mAh per pcs
Waktu pengisian	2 jam
Waktu penggunaan baterai	5 jam (maks. iluminasi) 20 jam (min. iluminasi)
Kontrol kecerahan	Tombol digital on/off dan kontrol kecerahan dengan layar LED
Magnifier optional	M250 2.5X (Hitam & Putih)/M350 3.5X (Hitam & Putih)/ M450 2.5-3.5X Hitam / M-4X, CM400/ DM500 /DM600
Kamera optional	1080p HD wireless camera