


# ENDO

## ENDO Operating Lamp with HD Camera



### Spesifikasi Operating Lamp

Tegangan Input	100 - 240 VAC
Konsumsi Daya	240 VA
Intensitas Cahaya	160.000 + 160.000 lux
Suhu Warna	3.500 – 5.500 K (dapat diatur)
Umur LED	>54.000
Fitur Tambahan	- Layar sentuh 5" TFT LCD - Mode Endo - Dapat diputar 360° - Handle dapat dilepas

### Spesifikasi HD Camera

Sensor	Exmor R CMOS 1/2.5" (8.51 MP)
Fitur	Image stabilization, freeze
Rotasi Kamera	360°