

ENDO

ENDO LED Operating Lamp EI.OL.8



Spesifikasi

Power supply	100-240VAC
Frekuensi	50-60Hz
Lamp head power	LED 4 lampu : 48W LED 5 lampu : 60W
Diameter lamp head	LED 4 lampu : 740mm LED 5 lampu : 770mm
Berat lamp head	LED 4 lampu : 15Kg LED 5 lampu : 17Kg
Pencahayaan pusat (Ec)	LED 4 lampu : 120klx LED 5 lampu : 160klx
Diameter bidang cahaya	220-320mm
Kedalaman bidang cahaya L1 + L2 (60%)	600mm
Color rendering index (CRI Ra)	Min. 95 (tipe 97) (R9: tipe 90)
Suhu warna	4000K
Penyesuaian cahaya elektrik (7 level)	LED 4 lampu : 45.000-120.000lx LED 5 lampu : 60.000-160.000lx
Masa pakai LED	≤60.000 jam