

ENDO

ENDO Electric Operating Table

EI.OT.3



Aksesoris

Spesifikasi Teknis

Tegangan Input	220 VAC \pm 10%
Frekuensi	50 Hz \pm 60 Hz
Beban Maksimal	250 kg

Dimensi Dasar

Panjang keseluruhan	2000 \pm 100mm
Lebar	500 \pm 20mm
Tinggi permukaan meja	750 \pm 50mm
Ketinggian pengangkatan meja	250 \pm 30mm
Kenaikan papan pinggang	100 \pm 20mm

Dimensi Dasar

Trendelenburg	$\geq 15^\circ$
Reverse Trendelenburg	$\geq 15^\circ$
Miring kiri	$\geq 15^\circ$
Miring kanan	$\geq 15^\circ$
Papan kepala naik	$\geq 40^\circ$
Papan kepala turun	$\geq 90^\circ$
Papan kaki turun	$\geq 90^\circ$
Papan kaki ke luar	$\geq 180^\circ$
Papan punggung naik	$\geq 45^\circ$
Papan punggung turun	$\geq 20^\circ$