

ENDO

ENDO Semi Automatic Urine Analyzer EI.UA



- 18 Parameter Leukocyte, Ketone, Nitrite, Urobilinogen, Bilirubin, Protein, Glucose, Specific Gravity, Blood, pH, Ascorbic Acid, Microalbumin, Creatinine, Calcium, ACR, PCR, Color, Turbidity
- Akurasi tinggi dan efisiensi tinggi
- Analyzer secara otomatis mendeteksi keberadaan strip, tipe dan memulai tes
- Baterai lithium terintegrasi

Spesifikasi	
Prinsip pengukuran	Reflectance photometry
Waktu reaksi	45 detik
Hasil	Hingga 550 T/jam
Display	Layar sentuh berwarna kapasitif resolusi tinggi 7 inci (1024*600)
Kapasitas memori	20.000 hasil
Interface	Port serial (Rs232), USB (pembaca barcode), wifi, bluetooth opsional
Power supply	Instrument: 12V DC; Adapter: 100V`240V, 50/60Hz
Dimensi	308mm x 210mm x 151mm (L x W x H)
Berat	4Kg