

ENDO Cannulated Drill EI.PT.CD



Parameter Teknis Utama



Tegangan	: 14 V	Pengisi Daya	
Daya	: $\geq 50W$	Daya keluaran	: AC 110-220V, 50/60Hz
Suhu	: $\leq 25^{\circ}C$	Tegangan Output	: 14V, 700M
Kecepatan(r.p.m)	: 0-350	Waktu pengisian	: ≤ 1.5 jam
		Daya Input	: 30VA
		Baterai	: Ni-MH 2/3A 1800mAh 9.6V